



# GENERAL SURGICAL INSTRUMENTS

Julespak International is a trusted Manufacturer & Exporter of Premium Surgical and Dental Instruments, proudly based in Sialkot, Pakistan — the world-renowned hub of precision surgical manufacturing.

With 18+ years of industry experience, we specialize in producing high-quality, durable, and internationally compliant medical instruments that meet the exact standards of hospitals, clinics, distributors, and healthcare brands worldwide.

Each product is crafted with precision engineering, premium-grade stainless steel, and strict quality control to ensure safety, performance, and long service life.

### **18 Years of Excellence & Trust**

**For nearly two decades, Julespak International has built a strong reputation for:**

**Consistent product quality**

**Reliable global exports**

**Custom manufacturing (OEM / Private Label)**

**Competitive pricing**

**On-time delivery**

Our experience enables us to understand international market requirements and deliver products that comply with global medical standards.

### **Why Choose Julespak International?**

- ✓ 18+ Years Manufacturing Expertise
- ✓ Factory-Direct Prices
- ✓ OEM & Custom Branding Available
- ✓ Skilled Workforce & Modern Machinery
- ✓ Strict Quality Inspection at Every Stage
- ✓ Worldwide Export Capability

### **Our Mission**

To provide world-class surgical and dental instruments that enhance medical performance, ensure patient safety, and build long-term partnerships with our global clients.

### **Our Vision**

To become a globally recognized name in surgical manufacturing, known for innovation, quality excellence, and customer satisfaction.



|                                                          |         |
|----------------------------------------------------------|---------|
| Bandage, Crown and Dissecting Scissors.....              | 04-11   |
| Episiotomy, Gynecological and Iris Scissors.....         | 12-14   |
| Angiotribes, Atrauma and Bile Duct Clamps.....           | 15-16   |
| Micro and Nasal Scissors.....                            | 17-20   |
| Plaster, Stitch and Wire Cutting Scissors.....           | 21-22   |
| SuperEdge SuperCut Scissors.....                         | 23-26   |
| Umbilical, Tonsil and Titanium Scissors.....             | 27-29   |
| Tungsten Carbide Scissors UltraCut.....                  | 30-33   |
| Tenotomy, Thorax and Vascular Scissors.....              | 34-36   |
| Operating Scissors.....                                  | 37-39   |
| Bulldog Clamps.....                                      | 40-42   |
| Dissecting and Ligature Forceps.....                     | 43-46   |
| Hemostatic Forceps.....                                  | 47-52   |
| Hysterectomy Forceps.....                                | 53-55   |
| Kidney Pedicle Clamps.....                               | 56      |
| Micro Vessel Clips and Nerve Approximators.....          | 57-58   |
| Amputation and Autopsy (Postmortem) Knives.....          | 59-61   |
| Skin Graft, Trigeminal, Tonsil and Cartilage Knives..... | 62-63   |
| Organ Grasping and Peritoneum Forceps.....               | 64-65   |
| Towel Clamps.....                                        | 66-68   |
| Tubing Clamps.....                                       | 69-71   |
| Dissecting and Operating Knives.....                     | 72-74   |
| Micro Scalpel and Scalpel Handles.....                   | 75-77   |
| Dressing Forceps.....                                    | 78-80   |
| Sponge Holding and Sinus Forceps.....                    | 81-83   |
| Atrauma, Cilia and Dissecting Forceps.....               | 84-87   |
| Diam-n-Dust Micro Ring Forceps.....                      | 88-90   |
| Diam-n-Dust Micro Suturing Forceps.....                  | 91      |
| Dressing and Micro Atrauma Forceps.....                  | 92-94   |
| Micro Ophthalmic Forceps.....                            | 95-98   |
| Micro Suturing and Splinter Forceps.....                 | 99-101  |
| UltraGripX Tungsten Carbide Forceps.....                 | 102-103 |
| Tissue Grasping and Tumor Grasping Forceps.....          | 104-106 |
| Percussion and Sensibility Instruments.....              | 107     |
| ENT Diagnostic Sets.....                                 | 108     |
| Standard Laryngoscopes.....                              | 109-111 |
| Fiber Optic Laryngoscopes.....                           | 112-114 |
| LED Fiber Optic Laryngoscopes.....                       | 115-117 |
| Delicate Hooks.....                                      | 118-120 |
| Lid Retractors & Saddle Hooks.....                       | 121-122 |
| Nasal Hooks.....                                         | 123     |
| Nerve Hooks.....                                         | 124-125 |
| Probes and Directors.....                                | 126-128 |
| Suction Tubes.....                                       | 129-132 |
| Trocars.....                                             | 133-134 |
| Wound Retractors.....                                    | 135-137 |
| General Retractors.....                                  | 138-150 |





Jl-32-101

**Smith-Mod. USA Bandage  
Scissor 18 cm - 7"**




Jl-32-102

**Smith-Mod. USA Bandage  
Scissor 23 cm - 9"**



Jl-32-103

**Smith-Mod. USA Bandage  
Scissor 20 cm - 8"**



Jl-32-104

**Bandage Scissor Straight 15 cm  
- 6"**



Jl-32-105

**Lorenz Bandage Scissor Heavy  
Pattern 23 cm - 9"**



Jl-32-106

**Lister Bandage Scissor 14.5 cm  
- 5 3/4"**



Jl-32-107

**Lister Bandage Scissor 18.5 cm  
- 7 1/4"**



Jl-32-108

**Lister Bandage Scissor 20 cm -  
8"**



Jl-32-109

**Lister Bandage Scissor One  
Large Ring 18.5 cm - 7 1/4"**





Jl-32-110

**Lister Bandage Scissor One  
Large Ring 21 cm - 8 1/4"**



Jl-32-111

**Lister-Excentric Bandage  
Scissor One Serrated Cutting  
Edge 16 cm - 6 1/4"**



Jl-32-112

**Lister-Excentric Bandage  
Scissor One Serrated Cutting  
Edge 19 cm - 7 1/2"**



Jl-32-113

**Esmarch Bandage Scissor  
Heavy Pattern - Short  
Blades 20 cm - 8"**



Jl-32-114

**Esmarch Bandage Scissor  
Heavy Pattern - Short  
Blades 23 cm - 9"**



Jl-32-115

**Universal Bandage Scissor  
Plastic Handle - Black 14.5 cm -  
5 3/4"**



Jl-32-116

**Universal Bandage Scissor  
Plastic Handle - Black 16 cm - 6  
1/4"**



Jl-32-117

**Universal Bandage Scissor  
Plastic Handle - Black 18 cm - 7  
1/2"**



Jl-32-118

**Universal Bandage Scissor  
Plastic Handle - Blue 14.5 cm -  
5 3/4"**



JI-32-119

**Universal Bandage Scissor**  
Plastic Handle - Blue 18 cm - 7  
1/2"



JI-32-120

**Universal Bandage Scissor**  
Plastic Handle - Blue 16 cm - 6  
1/4"



JI-32-121

**Universal Bandage Scissor**  
Plastic Handle - Red 14.5 cm - 5  
3/4"



JI-32-122

**Universal Bandage Scissor**  
Plastic Handle - Red 16 cm - 6  
1/4"




JI-32-123

**Universal Bandage Scissor**  
Plastic Handle - Red 18 cm - 7  
1/2"




JI-32-124

**Universal Bandage Scissor**  
Plastic Handle - Green 14.5 cm -  
5 3/4"




JI-32-125

**BeeBee Crown Scissor**  
Straight 10.5 cm - 4 1/4"



JI-32-126


**BeeBee Crown Scissor**  
Curved 10.5 cm - 4 1/4"



JI-32-127

**BeeBee Crown Scissor Straight -**  
One Toothed Cutting Edge 12  
cm - 4 3/4"





Jl-32-128

**BeeBee Crown Scissor Curved -  
One Toothed Cutting Edge 12  
cm - 4 3/4"**



Jl-32-129

**Mayo Dissecting Scissor  
Straight 14.5 cm - 5 3/4"**



Jl-32-130

**Mayo Dissecting Scissor  
Straight 15 cm - 6"**



Jl-32-131

**Mayo Dissecting Scissor  
Straight 17 cm - 6 3/4"**



Jl-32-132

**Mayo Dissecting Scissor  
Straight 23.5 cm - 9 1/4"**



Jl-32-133

**Mayo Dissecting Scissor Curved  
14.5 cm - 5 3/4"**



Jl-32-134

**Mayo Dissecting Scissor Curved  
15 cm - 6"**




Jl-32-135

**Mayo Dissecting Scissor Curved  
17 cm - 6 3/4"**



Jl-32-136

**Mayo Dissecting Scissor Curved  
23.5 cm - 9 1/4"**



Jl-32-137

**Mayo Dissecting Scissor  
Straight - With Chamfered  
Blades 19 cm - 7 1/2"**



Jl-32-138

**Mayo Dissecting Scissor  
Straight - With Chamfered  
Blades 21.5 cm - 8 1/2"**



Jl-32-139

**Mayo Dissecting Scissor Curved  
- With Chamfered Blades 19 cm  
- 7 1/2"**



Jl-32-140

**Mayo Dissecting Scissor Curved  
- With Chamfered Blades 21.5  
cm - 8 1/2"**



Jl-32-141

**Mayo-Harrington Dissecting  
Scissor Straight 23 cm - 9"**



Jl-32-142

**Mayo-Harrington Dissecting  
Scissor Curved 22.5 cm - 9"**



Jl-32-143

**Mayo-Harrington Dissecting  
Scissor Curved 28 cm - 11"**



Jl-32-144

**Mayo-Stillé Dissecting Scissor  
Straight 15 cm - 6"**



Jl-32-145

**Mayo-Stillé Dissecting Scissor  
Straight 17 cm - 6 3/4"**





Jl-32-146

**Mayo-Stille Dissecting Scissor**  
**Straight 19.5 cm - 7 3/4"**



Jl-32-147

**Mayo-Stille Dissecting Scissor**  
**Straight 21.5 cm - 8 1/2"**



Jl-32-148

**Mayo-Stille Dissecting Scissor**  
**Curved 15 cm - 6"**



Jl-32-149

**Mayo-Stille Dissecting Scissor**  
**Curved 17 cm - 6 3/4"**



Jl-32-150

**Mayo-Stille Dissecting Scissor**  
**Curved 19.5 cm - 7 3/4"**



Jl-32-151

**Mayo-Stille Dissecting Scissor**  
**Curved 21.5 cm - 8 1/2"**




Jl-32-152

**Mayo-Stille Dissecting Scissor**  
**Straight - With Chamfered**  
**Blades 14 cm - 5 1/2"**



Jl-32-153

**Mayo-Stille Dissecting Scissor**  
**Straight - With Chamfered**  
**Blades 17 cm - 6 3/4"**



Jl-32-154

**Mayo-Stille Dissecting Scissor**  
**Curved - With Chamfered**  
**Blades 14 cm - 5 1/2"**



Jl-32-155

**Mayo-Stille Dissecting Scissor**  
Curved - With Chamfered  
Blades 17 cm - 6 3/4"



Jl-32-156

**Mayo-Stille Dissecting Scissor**  
Straight 15 cm - 6"



Jl-32-157

**Mayo-Stille Dissecting Scissor**  
Curved 15 cm - 6"




Jl-32-158

**Mayo-Noble Gynecological**  
Scissor Straight 16.5 cm - 6  
1/2"



Jl-32-159

**Mayo-Noble Gynecological**  
Scissor Curved 16.5 cm - 6 1/2"




Jl-32-160

**Kelly Dissecting Scissor /**  
Operating Scissor Straight 16  
cm - 6 1/4"



Jl-32-161

**Kelly Dissecting Scissor /**  
Operating Scissor Curved 16 cm  
- 6 1/4"



Jl-32-162

**Joseph Dissecting Scissor /**  
Operating Scissor Straight 15  
cm - 6"



Jl-32-163

**Joseph Dissecting Scissor /**  
Operating Scissor Curved 15 cm  
- 6"





Jl-32-164

**Weller Dissecting Scissor Heavy**  
Pattern 28 cm - 11"



Jl-33-101

**Braun-Stadler Episiotomy  
Scissor 14.5 cm - 5 3/4"**



Jl-33-102

**Braun-Stadler Episiotomy  
Scissor 20 cm - 8"**



Jl-33-103

**Braun-Stadler Episiotomy  
Scissor 22 cm - 8 3/4"**



Jl-33-104

**Sims Gynecological Scissor  
Straight 20 cm - 8"**



Jl-33-105

**Sims Gynecological Scissor  
Straight 23 cm - 9"**



Jl-33-106

**Sims Gynecological Scissor  
Curved 20 cm - 8"**



Jl-33-107

**Sims Gynecological Scissor  
Curved 23 cm - 9"**



Jl-33-108

**Doyen Gynecological Scissor  
Straight 18.5 cm - 7 1/4"**



Jl-33-109

**Doyen Gynecological Scissor  
Curved 18.5 cm - 7 1/4"**





Jl-33-110

**Sims-Siebold Gynecological  
Scissor Straight 24.5 cm - 9  
3/4"**



Jl-33-111

**Sims-Siebold Gynecological  
Scissor Straight 24.5 cm - 9  
3/4"**



Jl-33-112

**Sims Gynecological Scissor  
Straight - Shar/Sharp 20 cm -  
8"**



Jl-33-113

**Sims Gynecological Scissor  
Curved - Blunt/Blunt 21 cm - 8  
1/4"**



Jl-33-114

**Wertheim Gynecological Scissor  
Curved 20 cm - 8"**



Jl-33-115

**Wertheim Gynecological Scissor  
Curved 23 cm - 9"**



Jl-33-116

**Iris Scissor Straight 9 cm - 3  
1/2"**



Jl-33-117

**Iris Scissor Straight 10 cm - 4"**



Jl-33-118

**Iris Scissor Curved 9 cm - 3 1/2"**



Jl-33-119

**Iris Scissor Curved 10 cm - 4"**



Jl-33-120

**Iris Scissor With Extra Flat  
Blades 10 cm - 4"**



Jl-33-121

**Bonn Iris Scissor Straight 9 cm  
- 3 1/2"**




Jl-33-122

**Bonn Iris Scissor Curved 9 cm -  
3 1/2"**



Jl-33-123

**Iris Scissor Straight 11.5 cm - 4  
1/2"**



Jl-33-124

**Iris Scissor Curved 11.5 cm - 4  
1/2"**



Jl-33-125

**Iris Scissor Curved to Side 11.5  
cm - 4 1/2"**



Jl-33-126

**Iris Scissor Angled 11.5 cm - 4  
1/2"**



Jl-33-127

**Iris Scissor Straight 12.5 cm -  
5"**



Jl-41-101

**Fergusson Angiotribe Curved 20 cm - 8"**



Jl-41-102

**Fergusson Angiotribe Curved 16 cm - 6 1/4"**



Jl-41-103

**Fergusson Angiotribe Straight 20 cm - 8"**



Jl-41-104

**Fergusson Angiotribe Straight 16 cm - 6 1/4"**



Jl-41-105

**Wertheim Atrauma Vaginal Clamp Right Angled 25 cm - 9 3/4"**



Jl-41-106

**Dick Atrauma Vaginal Clamp Angled 27.5 cm - 10 3/4"**



Jl-41-107

**Dick Atrauma Vaginal Clamp Angled 26.5 cm - 10 1/2"**



Jl-41-108

**Dick Atrauma Vaginal Clamp Angled 25.5 cm - 10"**



Jl-41-109

**Guyon Atrauma Kidney Pedicle Clamp Curved 24 cm - 9 1/2"**





JL-41-110

**Guyon Atrauma Kidney Pedicle  
Clamp Curved 23 cm - 9"**



JL-41-111

**Uro-Tangential Atrauma  
Tangential Forcep Fig. 2 27 cm  
- 10 3/4"**



JL-41-112

**Uro-Tangential Atrauma  
Tangential Forcep Fig. 1 25 cm  
- 9 3/4"**



JL-41-113

**Nissen Bile Duct Clamp Curved  
Downwards 21.5 cm - 8 1/2"**



JL-41-114

**Lahey (Sweet) Bile Duct Clamp  
Curved 23 cm - 9"**



JL-41-115

**Lahey Bile Duct Clamp  
Angled 22.5 cm - 8 3/4"**



JL-41-116

**Lahey (Sweet) Bile Duct Clamp  
Curved 19.5 cm - 7 3/4"**



Jl-34-101

**Jaboma Micro Scissor Angled -  
For Cleft Palate 18 cm - 7"  
Blade Size 12 mm**



Jl-34-102

**Jaboma Micro Scissor Angled -  
One Probed Tip 18 cm - 7"  
Blade Size 12 mm**



Jl-34-103

**Jaboma Micro Scissor Angled -  
One Probed Tip 18 cm - 7"  
Blade Size 12 mm**



Jl-34-104

**Micro Scissor Curved  
Downwards - Bayonet  
Shaped 18.5 cm - 7 1/4"**



Jl-34-105

**Micro Scissor Curved Upwards -  
Bayonet Shaped 18.5 cm - 7  
1/4"**



Jl-34-106

**Micro Scissor Straight - Bayonet  
Shaped 19 cm - 7 1/2"**



Jl-34-107

**Micro Scissor Curved - Bayonet  
Shaped 18.5 cm - 7 1/4"**



Jl-34-108

**Micro Scissor Straight - Bayonet  
Shaped 18.5 cm - 7 1/4"**



Jl-34-109

**Micro Scissor Curved - Bayonet  
Shaped 23 cm - 9"**



Jl-34-110

**Micro Scissor Straight - Bayonet  
Shaped 23 cm - 9"**



Jl-34-111

**Potts Micro Scissor Curved -  
Bayonet Shaped 23 cm - 9"**



Jl-34-112

**Jacobson Micro Scissor Straight  
- Bayonet Shaped 23 cm - 9"**



Jl-34-113

**Micro Scissor Curved - Flat  
Handle 23 cm - 9" Blade Size 10  
mm**



Jl-34-114

**Micro Scissor Curved - Flat  
Handle 21 cm - 8 1/4" Blade  
Size 10 mm**



Jl-34-115

**Micro Scissor Curved - Flat  
Handle 18 cm - 7" Blade Size 10  
mm**



Jl-34-116

**Micro Scissor Curved - Flat  
Handle 15 cm - 6" Blade Size 10  
mm**



Jl-34-117

**Micro Scissor Straight - Flat  
Handle 23 cm - 9" Blade Size 10  
mm**



Jl-34-118

**Micro Scissor Straight - Flat  
Handle 21 cm - 8 1/4" Blade  
Size 10 mm**





JI-34-119

**Micro Scissor Straight - Flat**  
Handle 18 cm - 7" Blade Size 10 mm



JI-34-120

**Micro Scissor Straight - Flat**  
Handle 15 cm - 6" Blade Size 10 mm



JI-34-121

**Micro Scissor Curved - Round**  
Handle 15 cm - 6" Blade Size 10 mm



JI-34-122

**Micro Scissor Curved - Round**  
Handle 18 cm - 7" Blade Size 10 mm



JI-34-123

**Micro Scissor Curved - Round**  
Handle 21 cm - 8 1/4" Blade Size 10 mm



JI-34-124

**Micro Scissor Curved - Round**  
Handle 23 cm - 9" Blade Size 10 mm



JI-34-125

**Cinelli Nasal Scissor**  
Curved 11.5 cm - 4 1/2"



JI-34-126

**Fomon Nasal Scissor Curved** 14 cm - 5 1/2"



JI-34-127

**Caplan Septum Scissor One**  
Toothed Cutting Edge 20 cm - 8"



Jl-34-128

**Fomon Tonsil Scissor 15 cm - 6"**



Jl-34-129

**Cottle Nasal Scissor 16 cm - 6 1/4"**



Jl-34-130

**Heymann Nasal Scissor 18.5 cm - 7 1/4"**



Jl-34-131

**Heymann Nasal Scissor One Toothed Cutting Edge 18.5 cm - 7 1/4"**



Jl-34-132

**Converse (Graefe) Nasal Scissor Angled - Blunt/Blunt 10.5 cm - 4 1/4"**



Jl-34-133

**Converse Nasal Scissor Angled - Sharp/Sharp 10.5 cm - 4 1/4"**



Jl-35-101

Brun's Plaster of Paris Shear One Serrated Cutting Edge 24 cm - 9 1/2"



Jl-35-102

Plaster of Paris Shear One Serrated Cutting Edge 20 cm - 8"



Jl-35-103

Seutin Plaster of Paris Shear 23 cm - 9"



Jl-35-104

Northbent Stitch Scissor 13 cm - 5"



Jl-35-105

Spencer Stitch Scissor 13 cm - 5"



Jl-35-106

Spencer Stitch Scissor 9 cm - 3 1/2"



Jl-35-107

Buck Stitch Scissor 14.5 cm - 5 3/4"



Jl-35-108

Buck Stitch Scissor 18 cm - 7"



Jl-35-109


Littauer Stitch Scissor 14 cm - 5 1/2"






Jl-35-110

**Eiselsberg Stitch Scissor 12 cm  
- 4 3/4"**



Jl-35-111

**Universal Wire Cutting Scissor  
Angled - One Toothed Cutting  
Edge 12 cm - 4 3/4"**



Jl-35-112

**Harvey Wire Cutting Scissor  
Straight - One Toothed Cutting  
Edge 12.5 cm - 5"**

Plaster, Stitch and Wire Cutting Scissors



JI-36-101

**SuperEdge™ Operating Scissor**  
Straight - Sharp/Blunt 13 cm -  
5"



JI-36-102

**SuperEdge™ Operating Scissor**  
Straight - Sharp/Blunt 14.5 cm -  
5 3/4"



JI-36-103

**SuperEdge™ Operating Scissor**  
Curved - Sharp/Blunt 13 cm -  
5"



JI-36-104

**SuperEdge™ Operating Scissor**  
Curved - Sharp/Blunt 14.5 cm -  
5 3/4"



JI-36-105

**Mayo SuperEdge™ Dissecting**  
Scissor Straight 14.5 cm - 5  
3/4"



JI-36-106

**Mayo SuperEdge™ Dissecting**  
Scissor Straight 15 cm - 6"



JI-36-107

**Mayo SuperEdge™ Dissecting**  
Scissor Straight 17 cm - 6 3/4"



JI-36-108

**Mayo SuperEdge™ Dissecting**  
Scissor Curved 14.5 cm - 5 3/4"



JI-36-109

**Mayo SuperEdge™ Dissecting**  
Scissor Curved 15 cm - 6"



Jl-36-110

**Mayo SuperEdge™ Dissecting Scissor Curved 17 cm - 6 3/4"**



Jl-36-111

**Mayo-Harrington Dissecting Scissor Straight 23 cm - 9"**



Jl-36-112

**Mayo-Harrington SuperEdge™ Dissecting Scissor Curved 22.5 cm - 9"**



Jl-36-113

**Mayo-Stille Dissecting Scissor Straight 15 cm - 6"**



Jl-36-114

**Mayo-Stille Dissecting Scissor Straight 17 cm - 6 3/4"**



Jl-36-115

**Mayo-Stille SuperEdge™ Dissecting Scissor Curved 15 cm - 6"**



Jl-36-116

**Mayo-Stille Dissecting Scissor Curved 17 cm - 6 3/4"**



Jl-36-117

**Sims SuperEdge™ Gynecological Scissor Straight 20 cm - 8"**



Jl-36-118

**Sims SuperEdge™ Gynecological Scissor Curved 20 cm - 8"**





Jl-36-119

**Sims SuperEdge™  
Gynecological Scissor  
Curved 23 cm - 9"**



Jl-36-120

**Sims-Siebold SuperEdge™  
Gynecological Scissor  
Straight 24.5 cm - 9 3/4"**



Jl-36-121

**Sims-Siebold Gynecological  
Scissor Curved 24.5 cm - 9  
3/4"**



Jl-36-122

**Kelly SuperEdge™ Dissecting  
Scissor / Opreating Scissor  
Straight 16 cm - 6 1/4"**



Jl-36-123

**Kelly SuperEdge™ Dissecting  
Scissor / Opreating Scissor  
Curved 16 cm - 6 1/4"**



Jl-36-124

**Joseph SuperEdge™ Dissecting  
Scissor / Opreating Scissor  
Straight 15 cm - 6"**



Jl-36-125

**Joseph SuperEdge™ Dissecting  
Scissor / Opreating Scissor  
Curved 15 cm - 6"**



Jl-36-126

**Metzenbaum SuperEdge™  
Dissecting Scissor / Opreating  
Scissor Straight -  
Blunt/Blunt 11.5 cm - 4 1/2"**



Jl-36-127

**Metzenbaum SuperEdge™  
Dissecting Scissor / Opreating  
Scissor Curved -  
Blunt/Blunt 11.5 cm - 4 1/2"**



Jl-36-128

**Metzenbaum Super Cut  
Dissecting Scissor / Opreating  
Scissor Straight -  
Blunt/Blunt 15.5 cm - 6"**



Jl-36-129

**Metzenbaum SuperEdge™  
Dissecting Scissor / Opreating  
Scissor Curved -  
Blunt/Blunt 15.5 cm - 6"**



Jl-36-130

**Metzenbaum SuperEdge™  
Dissecting Scissor Straight -  
Blunt/Blunt 14.5 cm - 5 3/4"**



Jl-36-131

**Metzenbaum SuperEdge™  
Dissecting Scissor Curved -  
Blunt/Blunt 14.5 cm - 5 3/4"**



Jl-36-132

**Metzenbaum-  
Nelson SuperEdge™ Dissecting  
Scissor Straight -  
Blunt/Blunt 18 cm - 7"**



Jl-36-133

**Metzenbaum-  
Nelson SuperEdge™ Dissecting  
Scissor Straight -  
Blunt/Blunt 20.5 cm - 8"**



Jl-36-134

**Metzenbaum-  
Nelson SuperEdge™ Dissecting  
Scissor Straight -  
Blunt/Blunt 23 cm - 9"**



Jl-36-135

**Metzenbaum-  
Nelson SuperEdge™ Dissecting  
Scissor Straight -  
Blunt/Blunt 26 cm - 10 1/4"**



Jl-36-136

**Metzenbaum-  
Nelson SuperEdge™ Dissecting  
Scissor Curved - Blunt/Blunt 18  
cm - 7"**



JI-37-101

**Modell USA Umbilical Cord  
Scissor 10.5 cm - 4 1/4"**



JI-37-102

**Schuhmacher Umbilical Cord  
Scissor 16 cm - 6 1/4"**



JI-37-103

**Busch Umbilical Cord Scissor 16  
cm - 6 1/4"**



JI-37-104

**Busch Umbilical Cord Scissor  
One Toothed Cutting Edge 16  
cm - 6 1/4"**



JI-37-105

**Prince Tonsil Scissor 18 cm -  
7"**



JI-37-106

**Dean Tonsil Scissor Angled  
Upwards 17.5 cm - 7"**



JI-37-107

**Dean Tonsil Scissor Toothed -  
Angled Upwards 17.5 cm - 7"**



JI-37-108

**Good Tonsil Scissor 19.5 cm - 7  
3/4"**



JI-37-109

**Micro Scissor Curved - Flat  
Handle 21 cm - 8 1/4" Blade  
Size 10 mm**





JI-37-110

**Micro Scissor Curved - Flat**  
Handle 18 cm - 7" Blade Size 10 mm



JI-37-111

**Micro Scissor Curved - Flat**  
Handle 15 cm - 6" Blade Size 10 mm



JI-37-112

**Micro Scissor Straight - Flat**  
Handle 21 cm - 8 1/4" Blade Size 10 mm



JI-37-113

**Micro Scissor Straight - Flat**  
Handle 18 cm - 7" Blade Size 10 mm



JI-37-115

**Micro Scissor Straight - Flat**  
Handle 15 cm - 6" Blade Size 10 mm



JI-37-116

**Micro Scissor Curved Round -**  
Handle 15 cm - 6" Blade Size 10 mm



JI-37-117

**Micro Scissor Curved Round -**  
Handle 18 cm - 7" Blade Size 10 mm



JI-37-118

**Micro Scissor Curved Round -**  
Handle 21 cm - 8 1/4" Blade Size 10 mm



JI-37-119

**Micro Scissor Straight - Round**  
Handle 15 cm - 6" Blade Size 10 mm



Jl-37-114

**Micro Scissor Straight - Round  
Handle 18 cm - 7" Blade Size 10  
mm**



Jl-37-120

**Micro Scissor Straight - Round  
Handle 21 cm - 8 1/4" Blade  
Size 10 mm**

Umbilical, Tonsil and Titanium Scissors



JI-38-101

**UltraCut™ TC Lister Bandage  
Scissor 18 cm - 7"**



JI-38-102

**UltraCut™ TC Locklin Gum  
Scissor Angled - S Shaped - One  
Toothed Cutting Edge 16 cm - 6  
1/4"**




JI-38-103

**UltraCut™ TC Locklin Gum  
Scissor Angled - One Toothed  
Cutting Edge 16 cm - 6 1/4"**




JI-38-104

**UltraCut™ TC Dean Gum Scissor  
S Shaped - One Toothed Cutting  
Edge 17.5 cm - 7"**




JI-38-105

**UltraCut™ TC Universal Wire  
Cutting Scissor Angled - One  
Toothed Cutting Edge 12.5 cm -  
5"**



JI-38-106

**UltraCut™ TC Kelly Operating  
Scissor Curved 16 cm - 6 1/4"**



JI-38-107

**UltraCut™ TC Kelly Operating  
Scissor Straight 16 cm - 6 1/4"**



JI-38-108


**UltraCut™ TC Ragnell Dissecting  
Scissor Curved 12.5 cm - 5"**



JI-38-109

**UltraCut™ TC Angled Scissor  
Angled 11 cm - 4 1/2"**





Jl-38-110

**UltraCut™ TC DeBakey Vascular  
Scissor Angled 60° 22 cm - 8  
3/4"**



Jl-38-111

**UltraCut™ TC DeBakey Vascular  
Scissor Angled 45° 23 cm - 9"**



Jl-38-112

**UltraCut™ TC DeBakey Vascular  
Scissor Angled 25° 23 cm - 9"**



Jl-38-113

**UltraCut™ TC DeBakey Vascular  
Scissor Angled 60° 18 cm - 7"**




Jl-38-114

**UltraCut™ TC DeBakey Vascular  
Scissor Angled 45° 18.5 cm - 7  
1/4"**



Jl-38-115

**UltraCut™ TC DeBakey Vascular  
Scissor Angled 25° 19 cm - 7  
1/2"**



Jl-38-116

**UltraCut™ TC Iris Scissor  
Curved 11.5 cm - 4 1/2"**



Jl-38-117

**UltraCut™ TC Iris Scissor  
Straight 11.5 cm - 4 1/2"**



Jl-38-118

**UltraCut™ TC Goldman-Fox Gum  
Scissor Curved 13 cm - 5"**



Jl-38-119

**UltraCut™ TC Goldman-Fox Gum  
Scissor Straight 13 cm - 5"**



Jl-38-120

**UltraCut™ TC Stevens Tenotomy  
Scissor Curved -  
Sharp/Sharp 11.5 cm - 4 1/2"**




Jl-38-121

**UltraCut™ TC Stevens Tenotomy  
Scissor Curved -  
Sharp/Sharp 10.5 cm - 4 1/4"**



Jl-38-122

**UltraCut™ TC Stevens Tenotomy  
Scissor Straight -  
Sharp/Sharp 10.5 cm - 4 1/4"**



Jl-38-123

**UltraCut™ TC Stevens Tenotomy  
Scissor Straight -  
Sharp/Sharp 11.5 cm - 4 1/2"**




Jl-38-124

**UltraCut™ TC Stevens Tenotomy  
Scissor Curved -  
Blunt/Blunt 11.5 cm - 4 1/2"**



Jl-38-125

**UltraCut™ TC Stevens Tenotomy  
Scissor Curved -  
Blunt/Blunt 10.5 cm - 4 1/4"**



Jl-38-126

**UltraCut™ TC Stevens Tenotomy  
Scissor Straight -  
Blunt/Blunt 11.5 cm - 4 1/2"**



Jl-38-127

**UltraCut™ TC Stevens Tenotomy  
Scissor Straight -  
Blunt/Blunt 10.5 cm - 4 1/4"**



Jl-38-128

**UltraCut™ TC**  
**Parametrium Hysterectomy**  
**Scissor Strongly Curved 22.5**  
**cm - 8 3/4"**



Jl-38-129

**UltraCut™ TC**  
**Parametrium Hysterectomy**  
**Scissor Curved 23 cm - 9"**




Jl-38-130

**UltraCut™ TC Sims Uterine**  
**Scissor Curved 23 cm - 9"**




Jl-38-131

**UltraCut™ TC Sims Uterine**  
**Scissor Straight 23 cm - 9"**



Jl-38-132

**UltraCut™ TC**  
**Wertheim Gynecological Scissor**  
**Curved 35 cm - 14"**




Jl-38-133

**UltraCut™ TC**  
**Wertheim Gynecological Scissor**  
**Curved 23 cm - 9"**



Jl-38-134

**UltraCut™ TC Freeman Face-lift**  
**Scissor Toothed 23 cm - 9"**



Jl-38-135


**UltraCut™ TC**  
**Freeman-Kaye Face-lift Scissor**  
**Toothed 18 cm - 7"**



Jl-38-136

**UltraCut™ TC Gorney Face-lift**  
**Scissor Toothed 18 cm - 7"**





JI-39-101

**Tenotomy Scissor Straight 11.5  
cm - 4 1/2"**



JI-39-102

**Tenotomy Scissor Curved 11.5  
cm - 4 1/2"**



JI-39-103

**Stevens Tenotomy Scissor  
Straight - Sharp/Sharp 10 cm -  
4"**



JI-39-104

**Stevens Tenotomy Scissor  
Straight - Sharp/Sharp 11 cm -  
4 1/2"**



JI-39-105

**Stevens Tenotomy Scissor  
Straight - Blunt/Blunt 10 cm -  
4"**



JI-39-106

**Stevens Tenotomy Scissor  
Straight - Blunt/Blunt 11 cm - 4  
1/2"**



JI-39-107

**Stevens Tenotomy Scissor  
Straight - Blunt/Blunt 12.5 cm -  
5"**



JI-39-108

**Stevens Tenotomy Scissor  
Curved - Sharp/Sharp 10 cm -  
4"**



JI-39-109

**Stevens Tenotomy Scissor  
Curved - Sharp/Sharp 11 cm - 4  
1/2"**



Jl-39-110

**Stevens Tenotomy Scissor**  
Curved - Blunt/Blunt 10 cm - 4"




Jl-39-111

**Stevens Tenotomy Scissor**  
Curved - Blunt/Blunt 11 cm - 4 1/2"




Jl-39-112

**Stevens Tenotomy Scissor**  
Curved - Blunt/Blunt 12.5 cm - 5"



Jl-39-113

**DeBakey Arteriotomy Scissor**  
Curved 17.5 cm - 7"



Jl-39-114

**Vascular Scissor Straight** 19 cm  
- 7 1/2"



Jl-39-115

**Vascular Scissor Curved** 19 cm -  
7 1/2"



Jl-39-116

**Diethrich Vascular Scissor**  
Angled 125° - With Guide 17.5  
cm - 7"



Jl-39-117

**DeBakey Vascular Scissor**  
Angled 45° 16.5 cm - 6 1/2"



Jl-39-118

**DeBakey Vascular Scissor**  
Angled 60° 16.5 cm - 6 1/2"



JI-39-119

**Potts-Smith Vascular Scissor**  
Angled 60° With Probe Tip 18  
cm - 7"



JI-39-120

**Potts-Smith Vascular Scissor**  
Angled 25° 19 cm - 7 1/2"



JI-39-121

**Potts-Smith Vascular Scissor**  
Angled 45° 18.5 cm - 7 1/4"




JI-39-122

**Potts-Smith Vascular Scissor**  
Angled 60° 18.5 cm - 7 1/4"



JI-39-123

**Potts-Smith Vascular Scissor**  
Angled on Flat 18.5 cm - 7  
1/4"



JI-39-124

**Diethrich-Potts Vascular Scissor**  
Angled 25° - Standard Blade 18  
cm - 7"





Jl-40-101

**Operating Scissor Straight -  
Blunt/Blunt 12 cm - 4 3/4"**



Jl-40-102

**Operating Scissor Straight -  
Blunt/Blunt 13 cm - 5"**




Jl-40-103

**Operating Scissor Straight -  
Blunt/Blunt 14.5 cm - 5 3/4"**



Jl-40-104

**Operating Scissor Straight -  
Blunt/Blunt 14.5 cm - 5 3/4"**



Jl-40-105

**Operating Scissor Straight -  
Blunt/Blunt 18.5 cm - 7 1/4"**



Jl-40-106

**Operating Scissor Straight -  
Blunt/Blunt 20.5 cm - 8"**



Jl-40-107

**Operating Scissor Curved -  
Blunt/Blunt 12 cm - 4 3/4"**




Jl-40-108

**Operating Scissor Curved -  
Blunt/Blunt 13 cm - 5"**



Jl-40-109

**Operating Scissor Curved -  
Blunt/Blunt 14.5 cm - 5 3/4"**



Jl-40-110

**Operating Scissor Curved -  
Blunt/Blunt 16.5 cm - 6 1/2"**



Jl-40-111

**Operating Scissor Curved -  
Blunt/Blunt 18.5 cm - 7 1/4"**



Jl-40-112

**Operating Scissor Curved -  
Blunt/Blunt 20.5 cm - 8"**



Jl-40-113

**Operating Scissor Straight -  
Sharp/Blunt 12 cm - 4 3/4"**



Jl-40-114

**Operating Scissor Straight -  
Sharp/Blunt 13 cm - 5"**



Jl-40-115

**Operating Scissor Straight -  
Sharp/Blunt 16.5 cm - 6 1/2"**



Jl-40-116

**Operating Scissor Straight -  
Sharp/Blunt 14.5 cm - 5 3/4"**



Jl-40-117

**Operating Scissor Straight -  
Sharp/Blunt 18.5 cm - 7 1/4"**



Jl-40-118

**Operating Scissor Straight -  
Sharp/Blunt 20.5 cm - 8"**



Jl-40-119

**Operating Scissor Curved -  
Sharp/Blunt 12 cm - 4 3/4"**



Jl-40-120

**Operating Scissor Curved -  
Sharp/Blunt 13 cm - 5"**



Jl-40-121

**Operating Scissor Curved -  
Sharp/Blunt 14.5 cm - 5 3/4"**



Jl-40-122

**Operating Scissor Curved -  
Sharp/Blunt 16.5 cm - 6 1/2"**



Jl-40-123

**Operating Scissor Curved -  
Sharp/Blunt 18.5 cm - 7 1/4"**



Jl-40-124

**Operating Scissor Curved -  
Sharp/Blunt 20.5 cm - 8"**





JI-42-101

**Dieffenbach Bulldog Clamp Straight**



JI-42-102

**Dieffenbach Bulldog Clamp Curved  
Stainless Steel, 6 cm - 2 1/4"**



JI-42-103

**Dieffenbach Bulldog Clamp Straight  
Stainless Steel, 3.5 cm - 1 1/2"**



JI-42-104

**Dieffenbach Bulldog Clamp Curved  
Stainless Steel, 3.5 cm - 1 1/2"**



JI-42-105

**Dieffenbach Bulldog Clamp Curved  
Stainless Steel, 5 cm - 2"**



JI-42-106

**Dieffenbach Bulldog Clamp Straight  
Stainless Steel, 5 cm - 2"**



JI-42-107

**John-Hopkins Bulldog Clamp 3.5  
cm - 1 1/2"**



JI-42-108

**John-Hopkins Bulldog Clamp 5  
cm - 2"**



JI-42-109

**John-Hopkins Bulldog Clamp 5.5  
cm - 2 1/8"**



Jl-42-110

**John-Hopkins Bulldog Clamp 6.5  
cm - 1/2"**



Jl-42-111

**John-Hopkins Bulldog Clamp 7  
cm - 2 3/4"**



Jl-42-112

**John-Hopkins Bulldog Clamp 7  
cm - 2 3/4"**



Jl-42-113

**John-Hopkins Bulldog Clamp 9  
cm - 3 1/2"**



Jl-42-114

**John-Hopkins Clip Applying  
Forcep 24 cm - 9 1/2"**



Jl-42-115

**Diethrich Bulldog Clamp  
Straight 45 mm Jaw Length 8  
mm**



Jl-42-116

**Diethrich Bulldog Clamp  
Straight 50 mm Jaw Length 12  
mm**



Jl-42-117

**Diethrich Bulldog Clamp  
Straight 60 mm Jaw Length 20  
mm**



Jl-42-118

**Diethrich Bulldog Clamp  
Angled 40 mm Jaw Length 8 mm**



Jl-42-119

**Diethrich Bulldog Clamp**  
**Angled 45 mm Jaw Length 12**  
**mm**



Jl-42-120

**Diethrich Bulldog Clamp**  
**Angled 50 mm Jaw Length 20**  
**mm**



Jl-42-121

**DeBakey-Diethrich Atrauma**  
**Bulldog Clamp Straight 45 mm**  
**Jaw Length 10 mm**



Jl-42-122

**DeBakey-Diethrich Atrauma**  
**Bulldog Clamp Straight 50 mm**  
**Jaw Length 14 mm**



Jl-42-123

**DeBakey-Diethrich Atrauma**  
**Bulldog Clamp Straight 60 mm**  
**Jaw Length 20 mm**



Jl-42-124

**DeBakey-Diethrich Atrauma**  
**Bulldog Clamp Straight 45 mm**  
**Jaw Length 10 mm**





Jl-43-101

**Gray Dissecting and Ligature  
Forcep Curved S Shaped 22 cm  
- 8 3/4"**



Jl-43-102

**Gray Dissecting and Ligature  
Forcep Curved S Shaped 22 cm  
- 8 3/4"**



Jl-43-103

**Finochietto Dissecting and  
Ligature Forcep Right Angled 24  
cm - 9 1/2"**



Jl-43-104

**Barre Dissecting and Ligature  
Forcep Delicate 28 cm - 11"**



Jl-43-105

**Barre Dissecting and Ligature  
Forcep Extra Delicate 28 cm -  
11"**



Jl-43-106

**Kantrowitz Dissecting and  
Ligature Forcep Right  
Angled 27.5 cm - 10 3/4"**



Jl-43-107

**Kantrowitz Dissecting and  
Ligature Forcep Right Angled 25  
cm - 9 3/4"**



Jl-43-108

**Kantrowitz Dissecting and  
Ligature Forcep Right Angled 20  
cm - 8"**



Jl-43-109

**Lahey Dissecting and Ligature  
Forcep Curved 18.5 cm - 7 1/4"**



Jl-43-110

**Wickstroem Dissecting and  
Ligature Forcep Right  
Angled 24.5 cm - 9 3/4"**



Jl-43-111

**Wickstroem Dissecting and  
Ligature Forcep Right Angled 21  
cm - 8 1/4"**



Jl-43-112

**Dissecting and Ligature Forcep  
Curved 26 cm - 10 1/4"**



Jl-43-113

**O-Shaugnessy Dissecting and  
Ligature Forcep Curved 23 cm -  
9"**



Jl-43-114

**O-Shaugnessy Dissecting and  
Ligature Forcep Curved 20 cm -  
8"**



Jl-43-115

**O-Shaugnessy Dissecting and  
Ligature Forcep Curved 18 cm -  
7"**



Jl-43-116

**O-Shaugnessy Dissecting and  
Ligature Forcep Curved 15 cm -  
6"**



Jl-43-117

**Mixer-  
O-Shaugnessy Dissecting and  
Ligature Forcep Curved 18.5 cm  
- 7"**



Jl-43-118

**Mc Quigg-Mixer Dissecting and  
Ligature Forcep Right Angled -  
Longitudinally Serrated 23.5 cm  
- 9 1/4"**



Jl-43-119

**Mixer Dissecting and Ligature  
Forcep Right Angled -  
Longitudinally Serrated 28 cm -  
11"**



Jl-43-120

**Mixer Dissecting and Ligature  
Forcep Right Angled -  
Longitudinally Serrated 26 cm -  
10 1/4"**



Jl-43-121

**Mixer Dissecting and Ligature  
Forcep Right Angled -  
Longitudinally Serrated 22.5 cm  
- 8 3/4"**



Jl-43-122

**Mixer Dissecting and Ligature  
Forcep Curved 23 cm - 9"**



Jl-43-123

**Mixer Dissecting and Ligature  
Forcep Curved 22.5 cm - 8 3/4"**



Jl-43-124

**Mixer Dissecting and Ligature  
Forcep Right Angled 18 cm - 7"**



Jl-43-125

**Mixer Dissecting and Ligature  
Forcep Right Angled 16 cm - 6  
1/4"**



Jl-43-126

**Mixer-Baby Dissecting and  
Ligature Forcep Curved 18.5 cm  
- 7 1/4"**



Jl-43-127

**Mixer-Baby Dissecting and  
Ligature Forcep Curved 13.5 cm  
- 5 1/4"**



JI-43-128

**Mixer-Baby Dissecting and Ligature Forcep Strongly Curved 13 cm - 5"**



JI-43-129

**Mixer-Baby Dissecting and Ligature Forcep Gently Curved 13.5 cm - 5 1/4"**



JI-43-130

**Zenker Dissecting and Ligature Forcep Curved 29.5 cm - 11 1/2"**



JI-43-131

**Zenker Dissecting and Ligature Forcep Curved 28.5 cm - 11 1/4"**



JI-43-132

**Kieback Dissecting and Ligature Forcep Curved 23.5 cm - 9 1/4"**



JI-43-133

**Rumel Dissecting and Ligature Forcep Curved 24 cm - 9 1/2"**



JI-43-134

**Rumel Dissecting and Ligature Forcep Curved 23.5 cm - 9 1/4"**



JI-43-135

**Rumel Dissecting and Ligature Forcep Curved 23 cm - 9"**



JI-43-136

**Overholt-Martin Dissecting and Ligature Forceps Fig. 1 18.5 cm - 7 1/4"**





Jl-44-101

**Pean Haemostatic Forcep 14.5 cm - 5 3/4"**



Jl-44-102

**Pean Haemostatic Forcep 16 cm - 6 1/4"**



Jl-44-103

**Kroenlein Haemostatic Forcep 13 cm - 5"**



Jl-44-104

**Kocher Haemostatic Forcep Straight - 1 x 2 Teeth 15 cm - 6"**



Jl-44-105

**Kocher Haemostatic Forceps Curved - 1 x 2 Teeth 15 cm - 6"**



Jl-44-106

**Spencer-Wells Haemostatic Forcep Straight 13 cm - 5"**



Jl-44-107

**Spencer-Wells Haemostatic Forcep Straight 14.5 cm - 5 3/4"**



Jl-44-108

**Spencer-Wells Haemostatic Forcep Straight 15.5 cm - 6"**



Jl-44-109

**Spencer-Wells Haemostatic Forcep Straight 18.5 cm - 7 1/4"**



Jl-44-110

**Spencer-Wells Haemostatic  
Forcep Straight 20 cm - 8"**



Jl-44-111

**Pean-Nippon Haemostatic  
Forceps Straight 14.5 cm - 5  
3/4"**



Jl-44-112

**Pean-Nippon Haemostatic  
Forceps Straight 16.5 cm - 6  
1/2"**



Jl-44-113

**Spencer-Wells Haemostatic  
Forcep Curved 13 cm - 5"**



Jl-44-114

**Spencer-Wells Haemostatic  
Forcep Curved 14 cm - 5 1/2"**



Jl-44-115

**Pean-Nippon Haemostatic  
Forceps Straight 18.5 cm - 7  
1/4"**



Jl-44-116

**Pean-Nippon Haemostatic  
Forceps Curved 14 cm - 5 1/2"**



Jl-44-117

**Spencer-Wells Haemostatic  
Forcep Curved 15 cm - 6"**



Jl-44-118

**Spencer-Wells Haemostatic  
Forcep Curved 18 cm - 7"**



JI-44-119

**Pean-Nippon Haemostatic  
Forceps Curved 16 cm - 6 1/4"**



JI-44-120

**Pean-Nippon Haemostatic  
Forceps Curved 18 cm - 7"**



JI-44-121

**Spencer-Wells Haemostatic  
Forcep Curved 20 cm - 8"**



JI-44-122

**Kocher-Nippon Haemostatic  
Forceps Straight - 1 x 2  
Teeth 14.5 cm - 5 3/4"**



JI-44-123

**Kocher-Nippon Haemostatic  
Forceps Straight - 1 x 2  
Teeth 16.5 cm - 6 1/2"**



JI-44-124

**Kocher-Nippon Haemostatic  
Forceps Straight - 1 x 2  
Teeth 18.5 cm - 7 1/4"**



JI-44-125

**Chaput Haemostatic Forcep 2 x  
3 Teeth 12.5 cm - 5"**



JI-44-126

**Tuffier Haemostatic Forcep 4 x  
5 Teeth 12.5 cm - 5"**



JI-44-127

**Kocher-Nippon Haemostatic  
Forceps Curved - 1 x 2  
Teeth 14.5 cm - 5 3/4"**



Jl-44-128

**Kocher-Nippon Haemostatic  
Forceps Curved - 1 x 2  
Teeth 16 cm - 6 1/4"**



Jl-44-129

**Kocher-Nippon Haemostatic  
Forceps Curved - 1 x 2  
Teeth 18 cm - 7"**



Jl-44-130

**Pean (Delicate) Haemostatic  
Forceps Straight 13.5 cm - 5  
1/4"**



Jl-44-131

**Pean (Delicate) Haemostatic  
Forceps Straight 14.5 cm - 5  
3/4"**



Jl-44-132

**Pean (Delicate) Haemostatic  
Forceps Straight 17 cm - 6  
3/4"**



Jl-44-133

**Pean (Delicate) Haemostatic  
Forceps Curved 13 cm - 5"**



Jl-44-134

**Pean (Delicate) Haemostatic  
Forceps Curved 14 cm - 5 1/2"**



Jl-44-135

**Pean (Delicate) Haemostatic  
Forceps Curved 16.5 cm - 6  
1/2"**



Jl-44-136


**Kocher (Delicate) Haemostatic  
Forceps Straight - 1 x 2  
Teeth 13 cm - 5"**






Jl-44-137

**Kocher (Delicate) Haemostatic Forceps Straight - 1 x 2 Teeth 14 cm - 5 1/2"**



Jl-44-138

**Pratt Haemostatic Forcep 15 cm - 6"**



Jl-44-139

**Martel-Rees Face-Lift Forcep 13.5 cm - 5 1/4"**



Jl-44-140

**Hartmann Haemostatic Forcep Straight 10 cm - 4"**



Jl-44-141

**Hartmann Haemostatic Forcep Curved 9.5 cm - 3 3/4"**



Jl-44-142

**Hartmann Haemostatic Forcep Straight - 1 x 2 Teeth 10 cm - 4"**



Jl-44-143

**Hartmann Haemostatic Forcep Curved - 1 x 2 Teeth 9.5 cm - 3 3/4"**



Jl-44-144

**Carrel Haemostatic Forcep 8.5 cm - 3 1/8"**



Jl-44-145

**Haemostatic Forcep Cross Serrated Jaw 12.5 cm - 5"**



JI-44-146

**Terrier Haemostatic Forcep 1 x  
2 Teeth 8.5 cm - 3 1/8"**



JI-44-147

**Halsted-Mosquito Haemostatic  
Forcep Straight 12.5 cm - 5"**



JI-44-148

**Halsted-Mosquito Haemostatic  
Forcep Straight 14.5 cm - 5  
3/4"**

Hemostatic Forceps



Jl-45-101

**Wertheim Hysterectomy Forcep**  
Curved - 1 x 2 Teeth 22.5 cm - 8  
3/4"



Jl-45-102

**Wertheim Hysterectomy Forcep**  
Curved - 1 x 2 Teeth 22.5 cm - 8  
3/4"



Jl-45-103

**Wertheim Hysterectomy Forcep**  
Curved 26 cm - 10 1/4"



Jl-45-104

**Moynihan Hysterectomy Forcep**  
Curved 23.5 cm - 9 1/4"



Jl-45-105

**Moynihan Hysterectomy Forcep**  
Curved 23 cm - 9"



Jl-45-106

**Burke Hysterectomy Forcep**  
Curved - Long Jaws 25 cm - 9  
3/4"



Jl-45-107

**Burke Hysterectomy Forcep**  
Curved - Short Jaws 25 cm - 9  
3/4"




Jl-45-108

**Kieback Hysterectomy Forcep**  
Curved 25 cm - 9 3/4"




Jl-45-109

**Berkeley-Bonney Hysterectomy**  
Forcep Curved 20 cm - 8"




Jl-45-110

**Kieback Hysterectomy Forcep**  
Straight 25 cm - 9 3/4"



Jl-45-111

**Gwilliam Hysterectomy Forcep**  
Curved 20 cm - 8"



Jl-45-112

**Gwilliam Hysterectomy Forcep**  
Straight 20 cm - 8"



Jl-45-113

**Holzbach Hysterectomy Forcep**  
Curved 24.5 cm - 9 3/4"



Jl-45-114

**Holzbach Hysterectomy Forcep**  
Straight 24.5 cm - 9 3/4"



Jl-45-115

**Kocher (Delicate) Haemostatic Forceps** Straight - 1 x 2  
Teeth 16 cm - 6 1/4"




Jl-45-116

**Second Hysterectomy Forcep**  
Curved 24 cm - 9 1/2"



Jl-45-117


**Second-Landau Hysterectomy Forcep** Straight 24 cm - 9 1/2"



Jl-45-118


**Wertheim Hysterectomy Forcep**  
Angled 22 cm - 8 3/4"






JI-45-119

**Wertheim**  
**(Rogers) Hysterectomy Forcep**  
Fig. 3 24.5 cm - 9 3/4"




JI-45-120

**Wertheim**  
**(Rogers) Hysterectomy Forcep**  
Fig. 2 25 cm - 9 3/4"




JI-45-121

**Wertheim**  
**(Rogers) Hysterectomy Forcep**  
Fig. 1 25.5 cm - 10"



JI-45-122

**Heaney (Rogers) Hysterectomy**  
**Forcep Fig. 4 20 cm - 8"**



JI-45-123

**Heaney (Rogers) Hysterectomy**  
**Forcep Fig. 3 20.5 cm - 8"**



JI-45-124

**Heaney (Rogers) Hysterectomy**  
**Forcep Fig. 2 21 cm - 8 1/4"**



Jl-46-101

**Herrick Kidney Pedicle Clamp**  
Double Angled - Longitudinally  
Serrated 23 cm - 9"



Jl-46-102

**Stille Kidney Pedicle Clamp**  
Curved - Grooved 23 cm - 9"



Jl-46-103

**Guyon Kidney Pedicle Clamp**  
Curved 24 cm - 9 1/2"



Jl-46-104

**Guyon Kidney Pedicle Clamp**  
Curved 21.5 cm - 8 1/2"



Jl-46-105

**Guyon Kidney Pedicle Clamp**  
Curved 19.5 cm - 7 3/4"



Jl-46-106

**Guyon Kidney Pedicle Clamp**  
Curved 24 cm - 9 1/2"



Jl-47-101

**Biemer Clip Applying Forcep  
Without Lock 14.5 cm - 5 3/4"**



Jl-47-102

**Biemer Clip Applying Forcep  
With Lock 14.5 cm - 5 3/4"**



Jl-47-103

**Biemer Vessel Clip 14 mm Jaw  
Opening - Jaw Length 4 mm - 4  
mm**



Jl-47-104

**Biemer Vessel Clip 17 mm Jaw  
Opening - Jaw Length 7 mm - 5  
mm**



Jl-47-105

**Biemer Vessel Clip 17 mm Jaw  
Opening - Jaw Length 7 mm - 5  
mm**




Jl-47-106

**Biemer-Muller Approximator  
Straight - Wide 17 mm Jaw  
Opening - Jaw Length 5 mm - 9  
mm**



Jl-47-107

**Biemer-Muller Approximator  
Angled - Wide 17 mm Jaw  
Opening - Jaw Length 5 mm - 9  
mm**




Jl-47-108

**Biemer-Muller Approximator  
Straight - Narrow 13 mm Jaw  
Opening - Jaw Length 3 mm -  
4.5 mm**



Jl-47-109

**Biemer-Muller Approximator  
Straight - Narrow 17 mm Jaw  
Opening - Jaw Length 5 mm - 9  
mm**



JI-47-110

**Biemer-Muller Approximator**  
Angled - Narrow 17 mm Jaw  
Opening - Jaw Length 5 mm - 9  
mm



JI-47-111

**Muller Approximator 17 mm**  
Diameter 1 mm



JI-47-112

**Muller Approximator 17 mm**  
Diameter 2 mm



JI-47-113

**Muller Approximator 17 mm**  
Diameter 3 mm



JI-47-114

**Muller Approximator 17 mm**  
Diameter 4 mm



JI-47-115

**Approximator Complete With**  
Key



JI-47-116

**Key For Approximator**



JI-47-117

**Approximator Only**



JI-47-118

**Nerve Approximator with**  
Needles Only





JI-51-101

**Liston Amputation Knife With  
Hollow Handle 26 cm - 10 1/4"  
Blade Size 130 mm**



JI-51-102

**Liston Amputation Knife With  
Hollow Handle 29 cm - 11 1/2"  
Blade Size 160 mm**



JI-51-103

**Liston Amputation Knife With  
Hollow Handle 32.5 cm - 12  
3/4" Blade Size 190 mm**



JI-51-104

**Liston Amputation Knife With  
Hollow Handle 35 cm - 13 3/4"  
Blade Size 220 mm**



JI-51-105

**Catlin Amputation Knife With  
Hollow Handle 29 cm - 11 1/2"  
Blade Size 160 mm**



JI-51-106

**Langenbeck Flap Knife 25 cm -  
9 3/4" Blade Size 120 mm**



JI-51-107

**Phalangeal Knife 23.5 cm - 9  
1/4" Blade Size 105 mm**



JI-51-108

**Interosseous Knife 24 cm - 9  
1/2" Blade Size 110 mm**



JI-51-109

**Resection Knife Solid Handle,  
Dissecting End 19 cm - 7 1/2"  
Blade Size 65mm**



Jl-51-110

**Langenbeck Resection Knife 18  
cm - 7" Blade Size 55mm**



Jl-51-111

**Langenbeck Resection Knife 18  
cm - 7" Blade Size 55mm**



Jl-51-112

**Virchow Autopsy Knife With  
Wooden Handle 25.5 cm - 10"  
Blade Size 100 mm**



Jl-51-113

**Virchow Cartilage Knife With  
Wooden Handle 21 cm - 8 1/4"  
Blade Size 80 mm**



Jl-51-114

**Walb Post Mortem Knife With  
Wooden Handle 26 cm - 10 1/4"  
Blade Size 110 mm**



Jl-51-115

**Walb Post Mortem Knife With  
Wooden Handle 29 cm - 11 1/2"  
Blade Size 140 mm**



Jl-51-116

**Walb Post Mortem Knife With  
Wooden Handle 31 cm - 12 1/4"  
Blade Size 170 mm**



Jl-51-117

**Virchow Brain Knife With Hollow  
Handle 26 cm - 10 1/4" Blade  
Size 160 mm**



Jl-51-118

**Virchow Brain Knife With Hollow  
Handle 33.5 cm - 13 1/4" Blade  
Size 200 mm**



JL-51-119

**Virchow Brain Knife With Hollow  
Handle 38 cm -15" Blade Size  
240 mm**

Amputation and Autopsy (Postmortem) Knives



Jl-54-101

**Humby Dermatome / Skin Graft  
Knife 32 cm - 13"**



Jl-54-102

**Watson Dermatome / Skin Graft  
Knife 30.5 cm - 12"**



Jl-54-103

**Dermatome Blade 15 cm - 6"**



Jl-54-104

**Schink Dermatome / Skin Graft  
Knife Adjustable From 0.1 to 2.0  
mm 31 cm - 12 1/4"**



Jl-54-105

**Schink Dermatome Spare  
Blade 10 cm - 4"**



Jl-54-106

**Schink Skin Fixing Plate**



Jl-54-107

**Silver Dermatome / Skin Graft  
Knife For Usual Razor Blade 19  
cm - 7 1/2"**



Jl-54-108

**Olivecrona Trigeminal  
Knife 20.5 cm - 8"**



Jl-54-109

**Trigeminal Knife 23 cm - 9"**





Jl-54-107

**Silver Dermatome / Skin Graft  
Knife For Usual Razor Blade 19  
cm - 7 1/2"**



Jl-54-108

**Olivecrona Trigeminal  
Knife 20.5 cm - 8"**



Jl-54-109

**Trigeminal Knife 23 cm - 9"**



Jl-54-110

**Fowler Tonsil Knife 20 cm - 8"**



Jl-54-111

**Abraham Tonsil Knife 21 cm - 8  
1/4"**



Jl-54-112

**Masing Cartilage Slitting  
Knife 13 cm - 5"**




JI-48-101

**Atrauma Organ Grasping /  
Haemorrhoidal Forcep 25 cm -  
9 3/4"**



JI-48-102

**Atrauma Organ Grasping /  
Haemorrhoidal Forcep 26 cm -  
10 1/4"**



JI-48-103

**Collin Organ Grasping Forcep 14  
cm - 5 1/2"**



JI-48-104

**Collin Tumor Grasping  
Forcep 17 cm - 6 3/4"**



JI-48-105

**Pennington Organ Grasping  
Forcep 15.5 cm - 6"**



JI-48-106

**Pennington Organ Grasping  
Forcep 21 cm - 8 1/4"**



JI-48-107

**Faure Peritoneum Forcep  
Curved - 1 x 2 Teeth 21 cm - 8  
1/4"**



JI-48-108

**Faure Peritoneum Forcep  
Curved - 1 x 2 Teeth 21 cm - 8  
1/4"**



JI-48-109

**Mikulicz Peritoneum Forcep  
Curved - 1 x 2 Teeth 20.5 cm -  
8"**



JI-48-110

**Mikulicz Peritoneum Forcep**  
Curved - Longitudinal  
Serrations - 1 x 2 Teeth 24.5 cm  
- 9 3/4"



JI-48-111

**Mikulicz Peritoneum Forcep**  
Curved - 1 x 2 Teeth 20 cm - 8"



JI-48-112

**Mikulicz Peritoneum Forcep**  
Curved - 1 x 2 Teeth 20 cm - 8"



JI-48-113

**Mikulicz Peritoneum Forcep**  
Curved - 1 x 2 Teeth 18.5 cm - 7  
1/4"



JI-48-114

**Mikulicz-Baby Peritoneum**  
Forcep Curved - 1 x 2  
Teeth 14.5 cm - 5 3/4"



Jl-49-101

**Lane Towel Clamp 19.5 cm - 7 3/4"**



Jl-49-102

**Lane Towel Clamp 15.5 cm - 6"**



Jl-49-103

**Lane Towel Clamp 13 cm - 5"**



Jl-49-104

**Towel Clamp For Paper Clothes 14 cm - 5 1/2"**



Jl-49-105

**Towel Clamp For Paper Clothes 11.5 cm - 4 1/2"**



Jl-49-106

**Peers-Bertram Towel Clamp With Jaw for Tube Fixation 14.5 cm - 5 3/4"**



Jl-49-107

**Non Perforating Towel Clamp 13.5 cm - 5 1/4"**



Jl-49-108

**Non Perforating Towel Clamp 10 cm - 4"**



Jl-49-109

**Tohoku Towel Clamp Blunt 13.5 cm - 5 1/4"**





Jl-49-110

**Tohoku Towel Clamp Blunt 10.5 cm - 4"**



Jl-49-111

**Bernhard Towel Clamp 16.5 cm - 6 1/2"**



Jl-49-112

**Lorna Towel Clamp Fine Serrated, Wide Jaw 13 cm - 5"**



Jl-49-113

**Lorna Towel Clamp Fine Serrated, Wide Jaw 9.5 cm - 3 3/4"**



Jl-49-114

**Roeder Towel Clamp 13.5 cm - 5 1/4"**



Jl-49-115

**Backhaus Towel Clamp 16 cm - 6 1/4"**



Jl-49-116

**Backhaus Towel Clamp 13 cm - 5"**



Jl-49-117

**Backhaus Towel Clamp 11.5 cm - 4 1/2"**



Jl-49-118

**Backhaus Towel Clamp 9 cm - 3 1/2"**



Jl-49-119

**Backhaus Towel Clamp With  
Tube Holder 13 cm - 5"**



Jl-49-120

**Robin Towel Clamp With Holders  
for Cable up to 13 mm 13.5 cm -  
5 1/4"**



Jl-49-121

**Robin Towel Clamp With Holders  
for Cable up to 5 mm 13.5 cm -  
5 1/4"**



Jl-49-122

**Jones-Martin Towel Clamp  
Sharp - With Cable Holder 8.5  
cm - 3 1/4"**



Jl-49-123

**Jones Towel Clamp Sharp 8.5  
cm - 3 1/4"**



Jl-49-124

**Jones Towel Clamp Sharp 5.5  
cm - 2 1/4"**



Jl-50-101

**Mohr Pinchcocks Width 80 mm**

Material: Stainless Steel  
Finish: Satin Finish (Mirror finish available upon request)



Jl-50-102

**Mohr Pinchcocks Width 70 mm**

Material: Stainless Steel  
Finish: Satin Finish (Mirror finish available upon request)



Jl-50-103

**Mohr Pinchcocks Width 60 mm**

Material: Stainless Steel  
Finish: Satin Finish (Mirror finish available upon request)



Jl-50-104

**Mohr Pinchcocks Width 50 mm**

Material: Stainless Steel  
Finish: Satin Finish (Mirror finish available upon request)



Jl-50-105

**Mohr Pinchcocks Width 40 mm**

Material: Stainless Steel  
Finish: Satin Finish (Mirror finish available upon request)



Jl-50-106

**Hoffmann Tubing Clamp Width 40 mm**

Material: Stainless Steel  
Finish: Satin Finish (Mirror finish available upon request)



Jl-50-107

**Hoffmann Tubing Clamp Width 30 mm**



Jl-50-108

**Hoffmann Tubing Clamp Width 25 mm**



Jl-50-109

**Hoffmann Tubing Clamp Width 25 mm**



Jl-50-110

**Hoffmann Tubing Clamp Width  
17 mm**



Jl-50-111

**Hoffmann Tubing Clamp Width  
13 mm**



Jl-50-112

**Brompton Tube Milking Roller**

Tubing Clamps



Jl-50-113

**Tubing Clamp With Safety  
Guard 20.5 cm - 8"**



Jl-50-114

**Tubing Clamp With Safety  
Guard 19 cm - 7 1/2"**



Jl-50-115

**Tubing Clamp With Safety  
Guard 15.5 cm - 6"**



Jl-50-116

**Tubing Clamp 20.5 cm - 8"**



Jl-50-117


**Tubing Clamp 19 cm - 7 1/2"**



Jl-50-118

**Tubing Clamp 15.5 cm - 6"**





JI-50-119

**Presbyterian Tubing Clamp 18  
cm - 7"**



JI-50-120

**Presbyterian Tubing Clamp 16.5  
cm - 6 1/2"**



JI-50-121

**Klinikum Berlin Tubing Clamp  
Extra Long Ratchet 14 cm - 5  
1/2"**

Tubing Clamps



JI-52-101

**Dissecting Knife / Operating Knife With Metal Handle 15 cm - 6" Blade Size 30 mm**



JI-52-102

**Dissecting Knife / Operating Knife With Metal Handle 15 cm - 6" Blade Size 36 mm**



JI-52-103

**Dissecting Knife / Operating Knife With Metal Handle 15.5 cm - 6" Blade Size 40 mm**



JI-52-104

**Dissecting Knife / Operating Knife With Metal Handle 16 cm - 6 1/4" Blade Size 45 mm**



JI-52-105

**Dissecting Knife / Operating Knife With Metal Handle 15 cm - 6" Blade Size 27 mm**



JI-52-106

**Dissecting Knife / Operating Knife With Metal Handle 15 cm - 6" Blade Size 35 mm**



JI-52-107

**Dissecting Knife / Operating Knife With Metal Handle 16 cm - 6 1/4" Blade Size 40 mm**



JI-52-108

**Dissecting Knife / Operating Knife With Metal Handle 16 cm - 6 1/4" Blade Size 47 mm**



JI-52-109

**Dissecting Knife / Operating Knife With Metal Handle 17 cm - 6 3/4" Blade Size 50 mm**



JI-52-110

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



JI-52-111

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



JI-52-112

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



JI-52-113

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



JI-52-114

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



JI-52-115

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



JI-52-116

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



JI-52-117

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



JI-52-118

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



Jl-52-119

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



Jl-52-120

**Dieffenbach Dissecting Knife /  
Operating Knife With Metal  
Handle 17 cm - 6 3/4"**



Jl-52-121

**Bergmann Dissecting Knife /  
Operating Knife With Metal  
Handle**



Jl-52-122

**Dissecting Knife / Operating  
Knife Bellied Blade - Fig. 1 14  
cm - 5 1/2"**



Jl-52-123

**Dissecting Knife / Operating  
Knife Bellied Blade - Fig. 3 14  
cm - 5 1/2"**



Jl-52-124

**Dissecting Knife / Operating  
Knife Bellied Blade - Fig. 4 14  
cm - 5 1/2"**





JI-53-101

**Scalpel Handle For Micro  
Scalpel Blades 13.5 cm - 5 1/4"**



JI-53-102

**Scalpel Handle For Micro  
Scalpel Blades 13.5 cm - 5 1/4"**



JI-53-103

**Caspar Scalpel Handle For Micro  
Scalpel Blades 18 cm - 7"**



JI-53-104

**Scalpel Handle For Micro  
Scalpel Blades 13 cm - 5"**



JI-53-105

**Landolt Scalpel Handle Bayonet  
Shaped, Round 25 cm - 10"**



JI-53-106

**Landolt Scalpel Handle Bayonet  
Shaped, Square Shank 21 cm -  
8 1/4"**



JI-53-107

**Bard-Parker Scalpel Handle No.  
3 Solid, 12 cm - 4 3/4"**



JI-53-108

**Bard-Parker Scalpel Handle No.  
4 Solid, 13.5 cm - 5 3/8 "**



JI-53-109

**Scalpel Handle No. 3 & 4 Solid,  
Double Ended 16 cm - 6 1/4"**



Jl-53-110

**Scalpel Handle No. 3 Solid,  
cm-Graduation 12 cm - 4 3/4".**



Jl-53-111

**Scalpel Handle No. 3L Solid,  
Long 21.5 cm - 8 1/2".**



Jl-53-112

**Scalpel Handle No. 4L Solid,  
Long 21 cm - 8 1/4".**



Jl-53-113

**Scalpel Handle No. 3L Solid,  
Angled, Long 21 cm - 8 1/4".**



Jl-53-114

**Scalpel Handle No. 4L Solid,  
Angled, Long 21 cm - 8 1/4".**



Jl-53-115

**Scalpel Handle No. 7 Solid 16  
cm - 6 1/4"**



Jl-53-116

**Scalpel Handle No. 7K Solid,  
Short 12.5 cm - 5"**



Jl-53-117

**Collin Scalpel Handle No. 8 With  
Hollow Handle 13.5 cm - 5 3/8 "**



Jl-53-118

**Collin Scalpel Handle No. 9 With  
Hollow Handle 14.5 cm - 5 3/4"**



Jl-53-119

**Scalpel Handle No. 3 With Round  
Hollow Handle 15.5 cm - 6"**



Jl-53-120

**Scalpel Handle No. 3 Solid,  
Round 14.5 cm - 5 3/4"**



Jl-53-121

**Scalpel Handle No. 3 Solid,  
Round, Angled 14 cm - 5 1/2"**

Micro Scalpel and Scalpel Handles



Jl-53-122

**Scalpel Handle No. 3 Solid,  
Round 15 cm - 6"**



Jl-55-101

**Ulrich Dressing Forceps  
Straight 22.5 cm - 8 3/4"**



Jl-55-102

**Ulrich Dressing Forceps  
Straight 26.5 cm - 10 1/2"**



Jl-55-103

**Ulrich Dressing Forceps Curved  
22.5 cm - 8 3/4"**



Jl-55-104

**Ulrich Dressing Forceps Curved  
26.5 cm - 10 1/2"**



Jl-55-105

**Bryant Dressing Forceps Curved  
13 cm - 5"**



Jl-55-106

**Martin Dressing Forceps  
Straight 17.5 cm - 7"**



Jl-55-107

**Martin Dressing Forceps Curved  
17 cm - 6 3/4"**



Jl-55-108

**Gross-Maier Dressing Forceps  
Straight - With Ratchet 20.5 cm  
- 8"**



Jl-55-109

**Gross-Maier Dressing Forceps  
Straight - With Ratchet 22 cm -  
8 3/4"**





Jl-55-110

**Gross-Maier Dressing Forceps  
Straight - With Ratchet 26.5 cm -  
10 1/2"**



Jl-55-111

**Gross-Maier Dressing Forceps  
Curved - With Ratchet 20 cm -  
8"**



Jl-55-112

**Gross-Maier Dressing Forceps  
Curved - With Ratchet 21.5 cm -  
8 1/2"**



Jl-55-113

**Gross-Maier Dressing Forceps  
Curved - With Ratchet 26 cm -  
10 1/4"**



Jl-55-114

**Gross-Maier Dressing Forceps  
Curved - Without Ratchet 20.5  
cm - 8"**



Jl-55-115

**Gross-Maier Dressing Forceps  
Straight - Without Ratchet 22  
cm - 8 3/4"**



Jl-55-116

**Gross-Maier Dressing Forceps  
Straight - Without Ratchet 26.5  
cm - 10 1/2"**



Jl-55-117

**Gross-Maier Dressing Forceps  
Curved - Without Ratchet 20 cm  
- 8"**



Jl-55-118

**Gross-Maier Dressing Forceps  
Curved - Without Ratchet 21.5  
cm - 8 1/2"**



Jl-55-119

**Gross-Maier Dressing Forcep**  
**Curved - Without Ratchet 26 cm**  
**- 10 1/4"**



Jl-55-120

**Pelkmann Foreign Body Forcep**  
**Straight 26 cm - 10 1/4"**



Jl-55-121

**Pelkmann Foreign Body Forcep**  
**Curved 26 cm - 10 1/4"**



Jl-55-122

**Tunneling Forcep Curved 39.5**  
**cm - 15 1/2"**



Jl-55-123

**Collin Dressing Forcep**  
**Straight 25.5 cm - 10"**



Jl-55-124

**Cheron Dressing Forcep**  
**Straight 23.5 cm - 9 1/4"**



JI-56-101

**Duplay Sponge Holding Forcep**  
**Straight 21 cm - 8 1/4"**



JI-56-102

**Duplay Sponge Holding Forcep**  
**Curved 21 cm - 8 1/4"**



JI-56-103

**Bozemann Sponge Holding**  
**Forcep Straight 26 cm - 10 1/4"**



JI-56-104

**Bozemann Sponge Holding**  
**Forcep Curved S Shaped 26 cm**  
**- 10 1/4"**



JI-56-105

**Bozemann-Douglas Sponge**  
**Holding Forcep Curved S**  
**Shaped - One Large Ring 26 cm**  
**- 10 1/4"**



JI-56-106

**Heywood-Smith Sponge Holding**  
**Forcep Straight 20 cm - 8"**



JI-56-107

**Heywood-Smith Sponge Holding**  
**Forcep Straight 21 cm - 8 1/4"**



JI-56-108

**Bergmann Sponge Holding**  
**Forcep Straight 24 cm - 9 1/2"**



JI-56-109

**Foerster Sponge Holding Forcep**  
**Straight - Smooth Jaw 25 cm - 9**  
**3/4"**



Jl-56-110

**Foerster Sponge Holding Forcep**  
Curved - Smooth Jaw 24.5 cm - 9 3/4"



Jl-56-111

**Foerster Sponge Holding Forcep**  
Straight 18 cm - 7"



Jl-56-112

**Foerster Sponge Holding Forcep**  
Straight 20.5 cm - 8"



Jl-56-113

**Foerster Sponge Holding Forcep**  
Straight 25 cm - 9 3/4"



Jl-56-114

**Foerster Sponge Holding Forcep**  
Curved 17.5 cm - 7"



Jl-56-115

**Foerster Sponge Holding Forcep**  
Curved 20 cm - 8"



Jl-56-116

**Foerster Sponge Holding Forcep**  
Curved 24.5 cm - 9 3/4"



Jl-56-117


**Foerster Sponge Holding Forcep**  
Straight 25 cm - 9 3/4"



Jl-56-118


**Foerster Sponge Holding Forcep**  
Curved 24.5 cm - 9 3/4"






JI-56-119

**Rampley Sponge Holding Forcep**  
Straight 18.5 cm - 7 1/4"



JI-56-120

**Rampley Sponge Holding Forcep**  
Straight 24.5 cm - 9 3/4"



JI-56-121

**Lister Sinus Forceps Straight** 13  
cm - 5"



JI-56-122

**Lister Sinus Forceps Straight** 16  
cm - 6 1/4"



JI-56-123

**Lister Sinus Forceps Straight** 18  
cm - 7"



Jl-57-101

**Cooley Atrauma Forcep 24 cm -  
9 1/2" Tip Size 2.0 mm**



Jl-57-102

**Cooley Atrauma Forcep 20 cm -  
8" Tip Size 2.0 mm**



Jl-57-103

**Cooley Atrauma Forcep 16 cm -  
6 1/4" Tip Size 2.0 mm**



Jl-57-104

**DeBakey Atrauma Forcep 30 cm -  
11 3/4" Tip Size 3.5 mm**



Jl-57-105

**DeBakey Atrauma Forcep 24.5  
cm - 9 3/4" Tip Size 3.5 mm**



Jl-57-106

**DeBakey Atrauma Forcep 20 cm -  
8" Tip Size 3.5 mm**



Jl-57-107

**DeBakey Atrauma Forcep 16 cm -  
6 1/4" Tip Size 3.5 mm**



Jl-57-108

**DeBakey Atrauma Forcep 30 cm -  
11 3/4" Tip Size 2.8 mm**



Jl-57-109

**DeBakey Atrauma Forcep 24.5  
cm - 9 3/4" Tip Size 2.8 mm**



Jl-57-110

**DeBakey Atrauma Forcep 20 cm**  
- 8" Tip Size 2.8 mm



Jl-57-111

**DeBakey Atrauma Forcep 16 cm**  
- 6 1/4" Tip Size 2.8 mm



Jl-57-112

**DeBakey Atrauma Forcep**  
Angled 29.5 cm - 11 1/2" Tip  
Size 2.0 mm



Jl-57-113

**Barraquer Cilia Forcep Smooth**  
Jaw 10.5 cm - 4"



Jl-57-114

**Beer Cilia Forcep Smooth Jaw 9**  
cm - 3 1/2"



Jl-57-115

**Douglas Cilia Forcep Wide Tip 9**  
cm - 3 1/2"



Jl-57-116

**Douglas Cilia Forcep Narrow**  
Tip 9 cm - 3 1/2"



Jl-57-117

**Bergh Cilia Forcep Serrated 8.5**  
cm - 3 1/4"



Jl-57-118

**Littauer Cilia Forcep Smooth**  
Jaw 8.5 cm - 3 1/4"



JI-57-119

**Standard Pattern Dissecting  
Forcep 1 x 2 Teeth 10.5 cm - 4  
1/4"**



JI-57-120

**Standard Pattern Dissecting  
Forceps 1 x 2 Teeth 11.5 cm - 4  
1/2"**



JI-57-121

**Standard Pattern Dissecting  
Forceps 1 x 2 Teeth 13 cm - 5"**



JI-57-122

**Standard Pattern Dissecting  
Forceps 1 x 2 Teeth 14.5 cm - 5  
3/4"**



JI-57-123

**Standard Pattern Dissecting  
Forcep 1 x 2 Teeth 16 cm - 6  
1/4"**



JI-57-124

**Standard Pattern Dissecting  
Forceps 1 x 2 Teeth 18 cm - 7"**



JI-57-125

**Standard Pattern Dissecting  
Forceps 1 x 2 Teeth 20 cm - 8"**



JI-57-126

**Standard Pattern Dissecting  
Forceps 1 x 2 Teeth 23 cm - 9"**



JI-57-127

**Standard Pattern Dissecting  
Forceps 1 x 2 Teeth 25 cm - 9  
3/4"**





JI-57-128

**Standard Pattern Dissecting  
Forceps 1 x 2 Teeth 30.5 cm -  
12"**



JI-57-129

**Dissecting Forceps 2 x 3  
Teeth 12 cm - 4 3/4"**



JI-57-130

**Dissecting Forceps 2 x 3  
Teeth 13 cm - 5"**

Atrauma, Cilia and Dissecting Forceps



Jl-58-101

**Diam-n-Dust™ Gerald Micro  
Ring Forcep Curved 18 cm - 7"  
Diameter 1.0 mm**



Jl-58-102

**Diam-n-Dust™ Gerald Micro  
Ring Forcep Straight 18 cm - 7"  
Diameter 1.0 mm**



Jl-58-103

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Short  
Platform 21 cm - 8 1/4"  
Diameter 1.0 mm Ø**



Jl-58-104

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Short  
Platform 18.5 cm - 7 1/4"  
Diameter 1.0 mm Ø**



Jl-58-105

**Diam-n-Dust™ Micro Ring  
Forcep Curved - With Long  
Platform 21 cm - 8 1/4"  
Diameter 1.0 mm Ø**



Jl-58-106

**Diam-n-Dust™ Micro Ring  
Forcep Curved - With Long  
Platform 18.5 cm - 7 1/4"  
Diameter 1.0 mm Ø**



Jl-58-107

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Long  
Platform 21 cm - 8 1/4"  
Diameter 1.0 mm Ø**



Jl-58-108

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Long  
Platform 18.5 cm - 7 1/4"  
Diameter 1.0 mm Ø**



Jl-58-109

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Counter  
Balance 18.5 cm - 7 1/4"  
Diameter 2.0 mm Ø**



JI-58-110

**Diam-n-Dust™ Micro Ring  
Forcep Curved - With Counter  
Balance 25 cm - 9 3/4"  
Diameter 1.0 mm Ø**



JI-58-111

**Diam-n-Dust™ Micro Ring  
Forcep Curved - With Counter  
Balance 23 cm - 9" Diameter  
1.0 mm Ø**



JI-58-112

**Diam-n-Dust™ Micro Ring  
Forcep Curved - With Counter  
Balance 21 cm - 8 1/4"  
Diameter 1.0 mm Ø**



JI-59-101

**Diam-n-Dust™ Micro Ring  
Forcep Curved - With Counter  
Balance 18.5 cm - 7 1/4"  
Diameter 1.0 mm Ø**



JI-59-102

**Diam-n-Dust™ Micro Ring  
Forcep Curved - With Counter  
Balance 15 cm - 6" Diameter  
1.0 mm Ø**



JI-59-103

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Counter  
Balance 25 cm - 9 3/4"  
Diameter 1.0 mm Ø**



JI-59-104

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Counter  
Balance 23 cm - 9" Diameter  
1.0 mm Ø**



JI-59-105

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Counter  
Balance 21 cm - 8 1/4"  
Diameter 1.0 mm Ø**



JI-59-106

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Counter  
Balance 18.5 cm - 7 1/4"  
Diameter 1.0 mm Ø**



Jl-59-107

**Diam-n-Dust™ Micro Ring  
Forcep Straight - With Counter  
Balance 15 cm - 6" Diameter  
1.0 mm Ø**



Jl-59-108

**Diam-n-Dust™ Micro Ring  
Forcep 15 cm - 6" Diameter 1 x  
0.5 mm Ø**



Jl-59-109

**Diam-n-Dust™ Micro Ring  
Forcep 21 cm - 8 1/4" Tip Size  
2.0 x 1.0 mm**

Diam-n-Dust Micro Ring Forceps



Jl-59-110

**Diam-n-Dust™ Micro Ring  
Forcep 18 cm - 7" Tip Size 2.0 x  
1.0 mm**



Jl-59-111

**Diam-n-Dust™ Micro Ring  
Forcep 15 cm - 6" Tip Size 2.0 x  
1.0 mm**



Jl-59-112

**Diam-n-Dust™ Micro Ring  
Forcep 23 cm - 9" Tip Size 2.0 x  
1.0 mm Ø**





Jl-60-101



Jl-60-102



Jl-60-103



Jl-60-104



Jl-60-105



Jl-60-106



Jl-60-107



Jl-60-108



Jl-60-109



Jl-60-110



Jl-60-111



Jl-60-112



JI-61-101

**Standard Pattern Dressing  
Forceps 10.5 cm - 4 1/4"**



JI-61-102

**Standard Pattern Dressing  
Forceps 11.5 cm - 4 1/2"**



JI-61-103

**Standard Pattern Dressing  
Forceps 13 cm - 5"**



JI-61-104

**Standard Pattern Dressing  
Forceps 14.5 cm - 5 3/4"**



JI-61-105

**Standard Pattern Dressing  
Forceps 16 cm - 6 1/4"**



JI-61-106

**Standard Pattern Dressing  
Forceps 20 cm - 8"**



JI-61-107

**Standard Pattern Dressing  
Forceps 23 cm - 9"**



JI-61-108

**Standard Pattern Dressing  
Forceps 25 cm - 9 3/4"**



JI-61-109

**Standard Pattern Dressing  
Forceps 30.5 cm - 12"**



Jl-61-110

**Slender Pattern Dressing  
Forceps 12 cm - 4 3/4"**



Jl-61-111

**Slender Pattern Dressing  
Forceps 13 cm - 5"**



Jl-61-112

**Slender Pattern Dressing  
Forceps 14.5 cm - 5 3/4"**



Jl-61-113

**Gerald (DeBakey) Micro  
Atrauma Forcep 20 cm - 8" Tip  
Size 1.0 mm**



Jl-61-114

**Gerald (DeBakey) Micro  
Atrauma Forcep 18 cm - 7" Tip  
Size 1.0 mm**



Jl-61-115

**Micro Atrauma Forcep 21 cm - 8  
1/4" Tip Size 1.2 mm**



Jl-61-116

**Micro Atrauma Forcep 18.5 cm -  
7 1/4" Tip Size 1.2 mm**



Jl-61-117

**Micro Atrauma Forcep With  
Counter Balance 21 cm - 8 1/4"  
Tip Size 1.2 mm**



Jl-61-118

**Micro Atrauma Forcep With  
Counter Balance 18.5 cm - 7  
1/4" Tip Size 1.2 mm**



Jl-61-119

**Micro Atrauma Forcep 21 cm - 8  
1/4" Tip Size 1.2 mm**



Jl-61-120

**Micro Atrauma Forcep 18 cm -  
7" Tip Size 1.2 mm**



Jl-61-121

**Micro Atrauma Forcep 15 cm -  
6" Tip Size 1.2 mm**



Jl-61-122

**Micro Atrauma Forcep 21 cm - 8  
1/4" Tip Size 1.2 mm**



Jl-61-123

**Micro Atrauma Forcep 18 cm -  
7" Tip Size 1.2 mm**



Jl-61-124

**Micro Atrauma Forcep 15 cm -  
6" Tip Size 1.2 mm**





Jl-62-101

**Desmarres Entropion Forcep  
Left 10 cm - 4"**



Jl-62-102

**Desmarres Entropion Forcep  
Right 10 cm - 4"**



Jl-62-103

**Wies Chalazion Forcep Grooved  
Upper Jaw 9.5 cm - 3 3/4"**



Jl-62-104

**Desmarres Chalazion Forcep  
Large - Oval Jaws 9 cm - 3 1/2"  
Jaw Size 20 x 13 mm**



Jl-62-105

**Francis Chalazion Forcep  
Medium - Oval Jaws 9.5 cm - 3  
3/4" Jaw Size 14 x 11 mm**



Jl-62-106

**Baird Chalazion Forcep Small -  
Oval Jaws 9.5 cm - 3 3/4" Jaw  
Size 11 x 8 mm**



Jl-62-107

**Hunt Chalazion Forcep Round  
Jaws 9.5 cm - 3 3/4" Jaw  
Diameter 12 mm Ø**



Jl-62-108

**Lambert Chalazion Forcep  
Round Jaws 9 cm - 3 1/2" Jaw  
Diameter 8 mm Ø**



Jl-62-109

**Downes Lid Clamp 9.5 cm - 3  
3/4"**



Jl-62-110

**Erhardt Lid Forcep 10.5 cm - 4"**



Jl-62-111

**Berke Ptosis Forcep  
Longitudinally Grooved Jaws -  
With Side Lock 11 cm - 4 1/4"  
Jaws Length 27 mm**



Jl-62-112

**Berke Ptosis Forcep  
Longitudinally Grooved Jaws -  
With Side Lock 11 cm - 4 1/4"  
Jaws Length 20 mm**



Jl-62-113

**Jameson Muscle Forcep Left -  
With Slide Lock 6 Teeth - Adult  
Size 10 cm - 4"**



Jl-62-114

**Jameson Muscle Forcep Right -  
With Slide Lock 6 Teeth - Adult  
Size 10 cm - 4"**



Jl-62-115

**Jameson Muscle Forcep Left -  
With Slide Lock 4 Teeth - Child  
Size 9.5 cm - 3 3/4"**



Jl-62-116

**Jameson Muscle Forcep Right -  
With Slide Lock 4 Teeth - Child  
Size 9.5 cm - 3 3/4"**



Jl-62-117

**Clayman Lens Holding Forcep  
Delicate Angled Jaws - Gently  
Curved 12.5 cm - 5"**



Jl-62-118

**Shepard Lens Holding Forcep  
Gently Curved - Duckbill Shaped  
Jaws 12.5 cm - 5"**



Jl-62-119

**Blaydes Lens Holding Forcep**  
Angled 45° - Very Delicate  
Narrow Jaws 11.5 cm - 4 1/2"



Jl-62-120

**Kratz Lens Holding Forcep Very**  
Delicate Narrow Jaws 12.5 cm -  
5"



Jl-62-121

**Bechert Lens Holding Forcep**  
Very Delicate Flat Jaws 12.5 cm  
- 5"



Jl-62-122

**Buratto Soft IOL Inserting**  
Forcep Extra Delicate Highly  
Polished Convex Jaws 10 cm -  
4"



Jl-62-123

**Faulkner Lens Holding Forcep**  
Round - Smooth Jaws 10 cm - 4"



Jl-62-124

**Livernois-McDonald IOL**  
Inserting Forcep Concave  
Jaws 10.5 cm - 4" Jaws Length 5  
mm



Jl-62-125

**Livernois IOL Holding Forcep**  
Curved - Highly Polished Round  
Jaws for Holding Soft IOLs 10  
cm - 4"



Jl-62-126

**Livernois IOL Holding Forcep**  
Straight - Highly Polished Round  
Jaws for Holding Soft IOLs 10.5  
cm - 4"



Jl-62-127

**Steinert IOL Inserting Forcep**  
Conically Shaped Jaws 10 cm -  
4"



JI-62-128

**Livernois IOL Folding Forcep**  
Concave Jaw Surfaces - For Soft  
IOLs 10.5 cm - 4"



JI-62-129

**Dodick IOL Folding Forcep For**  
Soft IOLs 11 cm - 4 1/4"



JI-62-130

**Rosenwasser Donor Lamella**  
Inserting Forcep With Vaulted  
Shanks 11 cm - 4 1/4"



JI-62-131

**Alfonso Nucleus Grasping**  
Forcep With Two Short Rows of  
Delicate Interlocking Teeth 10.5  
cm - 4"



JI-62-132

**Dodick Nucleus Cracker Cross**  
Action - Longitudinally Serrated  
Paddle Shaped Jaws 11 cm - 4  
1/4"



JI-62-133

**Ernest Nucleus Cracker Extra**  
Delicate Cross Action Serrated  
Paddle Tips 11 cm - 4 1/4"



JI-62-134

**Akahosii Piraco PreChopper**  
Delicate Jaws with Sharp  
Tips 11 cm - 4 1/4" Jaws Length  
4.5 mm



JI-62-135

**Arruga Capsule Forcep** 11 cm -  
4 1/4"



JI-62-136

**Jaffe Utility Forcep Curved -**  
Very Delicate Smooth Jaws 11  
cm - 4 1/4"





JI-63-101

**Gerl Suture Tying Forcep**  
Smooth Jaws and Curved  
Shanks 8.5 cm - 3 1/4" Jaws  
Length 6 mm



JI-63-102

**Gerl Suture Tying Forcep**  
Smooth Jaws and Curved  
Shanks 8.5 cm - 3 1/4" Jaws  
Length 6 mm



JI-63-103

**Fine Suture Tying Forcep**  
Curved 10.5 cm - 4"



JI-63-104

**Helveston Suture Tying Forcep**  
Straight 10.5 cm - 4"



JI-63-105

**Harms Suture Tying Forcep**  
Curved 10.5 cm - 4"



JI-63-106

**Harms Suture Tying Forcep**  
Straight 10.5 cm - 4"



JI-63-107

**Tennant Suture Tying Forcep**  
Curved - Round Handle with  
Guide Pin - Extra Delicate  
Smooth Jaws 10.5 cm - 4" Jaws  
Length 6 mm



JI-63-108

**Tennant Suture Tying Forcep**  
Curved - Round Handle with  
Guide Pin - Extra Delicate  
Smooth Jaws 10.5 cm - 4" Jaws  
Length 6 mm



JI-63-109

**Tennant Suture Tying Forcep**  
Straight - Round Handle with  
Guide Pin - Extra Delicate  
Smooth Jaws 10.5 cm - 4" Jaws  
Length 6 mm



Jl-63-110

**Tennant Suture Tying Forcep**  
Straight - Round Handle with  
Guide Pin - Extra Delicate  
Smooth Jaws 10.5 cm - 4" Jaws  
Length 6 mm



Jl-63-111

**Jaffe Suture Tying Forcep**  
Curved - Very Delicate Smooth  
Jaws 10.5 cm - 4" Jaws Length 6  
mm




Jl-63-112

**Jaffe Suture Tying Forcep**  
Straight - Very Delicate Smooth  
Jaws 10.5 cm - 4" Jaws Length 6  
mm



Jl-63-113

**Ralk Splinter Forcep** 15 cm - 6"




Jl-63-114

**Ralk Splinter Forcep** 15 cm - 6"



Jl-63-115

**Peet Splinter Forcep** 12.5 cm -  
5"



Jl-63-116

**Stieglitz Splinter Forcep** 13.5  
cm - 5 1/4"



Jl-63-117

**Stieglitz Splinter Forcep** 14 cm -  
5 1/2"



Jl-63-118

**Shaaf Splinter Forcep** 9.5 cm - 3  
3/4"



Jl-63-119

**Hunter Splinter Forcep 10.5 cm  
- 4"**



Jl-63-120

**Hunter Splinter Forcep 10.5 cm  
- 4"**



Jl-63-121

**Feilchenfeld Splinter  
Forcep 11.5 cm - 4 1/2"**



Jl-63-122

**Feilchenfeld Splinter Forcep 9  
cm - 3 1/2"**



Jl-63-123

**Feilchenfeld Splinter Forcep 7.5  
cm - 3 1/4"**



Jl-63-124

**Splinter Forcep Curved -  
Serrated Jaws 14.5 cm - 5 3/4"**



JI-65-101

**UltraGrip<sup>™</sup> TC**  
**DeBakey Atrauma Forcep 24.5**  
**cm - 9 3/4" Tip Size 2.0 mm**



JI-65-102

**UltraGrip<sup>™</sup> TC**  
**DeBakey Atrauma Forcep 20 cm**  
**- 8" Tip Size 2.0 mm**



JI-65-103

**UltraGrip<sup>™</sup> TC**  
**DeBakey Atrauma Forcep 16 cm**  
**- 6 1/4" Tip Size 2.0 mm**



JI-65-104

**UltraGrip<sup>™</sup> TC Micro-**  
**DeBakey Atrauma Forcep 24.5**  
**cm - 9 3/4" Tip Size 1.6 mm**



JI-65-105

**UltraGrip<sup>™</sup> TC Micro-**  
**DeBakey Atrauma Forcep 20 cm**  
**- 8" Tip Size 1.6 mm**



JI-65-106

**UltraGrip<sup>™</sup> TC Micro-**  
**DeBakey Atrauma Forcep 16 cm**  
**- 6 1/4" Tip Size 1.6 mm**



JI-65-107

**UltraGrip<sup>™</sup> TC**  
**Wangenstein Dressing**  
**Forcep 23 cm - 9"**



JI-65-108

**UltraGrip<sup>™</sup> TC**  
**Wangenstein Dressing**  
**Forcep 15 cm - 6"**



JI-65-109

**UltraGrip<sup>™</sup> TC**  
**Ferris-Smith Dissecting Forcep**  
**2 x 3 Teeth 17.5 cm - 7"**





Jl-65-110

**UltraGrip™ TC Dehler Dissecting  
Forcep 1 x 2 Teeth 25 cm - 9  
3/4"**



Jl-65-111

**UltraGrip™ TC Dehler Dissecting  
Forcep 1 x 2 Teeth 20 cm - 8"**



Jl-65-112

**UltraGrip™ TC Dehler Dissecting  
Forcep 1 x 2 Teeth 18 cm - 7"**



Jl-65-113

**UltraGrip™ TC Dehler Dissecting  
Forcep 1 x 2 Teeth 16 cm - 6  
1/4"**



Jl-65-114

**UltraGrip™ TC Dehler Dissecting  
Forcep 1 x 2 Teeth 14.5 cm - 5  
3/4"**



Jl-65-115

**UltraGrip™ TC  
Potts-Smith Dressing Forcep 25  
cm - 9 3/4"**

UltraGrip™ Tungsten Carbide Forceps



JI-64-101

**Selman Dissecting Forceps 15.5 cm - 6"**



JI-64-102

**Russian Modell Dissecting Forceps 15 cm - 6"**



JI-64-103

**Russian Modell Dissecting Forceps 20 cm - 8"**



JI-64-104

**Russian Modell Dissecting Forceps 25.5 cm - 10"**



JI-64-105

**Durante Dissecting Forceps 14.5 cm - 5 3/4"**



JI-64-106

**Durante Dissecting Forceps 16 cm - 6 1/4"**



JI-64-107

**Durante Dissecting Forceps 20.5 cm - 8"**



JI-64-108

**Stille-Barraya Dissecting Forceps 18 cm - 7"**



JI-64-109

**Stille-Barraya Dissecting Forceps 20 cm - 8"**



Jl-64-110

**Stille-Barraya Dissecting  
Forceps 25 cm - 9 3/4"**



Jl-64-111

**Mayo (Russian) Dissecting  
Forceps 23 cm - 9"**



Jl-64-112

**Singley-Tuttle Dissecting  
Forceps 18 cm - 7"**



Jl-64-113

**Samii Tumor Grasping Forcep  
Bayonet Shaped 24 cm - 9 1/2"  
Tip Size 3.0 mm Ø**



Jl-64-114

**Samii Tumor Grasping Forcep  
Bayonet Shaped 24 cm - 9 1/2"  
Tip Size 3.0 mm Ø**



Jl-64-115

**Yasargil Tumor Grasping Forcep  
Bayonet Shaped 22 cm - 8 3/4"  
Tip Size 3.0 mm Ø**



Jl-64-116

**Yasargil Tumor Grasping Forcep  
Bayonet Shaped 22 cm - 8 3/4"  
Tip Size 5.0 mm Ø**



Jl-64-117

**Yasargil Tumor Grasping  
Forceps Bayonet Shaped 22 cm  
- 8 3/4" Tip Size 3.0 mm Ø**



Jl-64-118

**Yasargil Tumor Grasping  
Forceps Bayonet Shaped 22 cm  
- 8 3/4" Tip Size 5.0 mm Ø**



Jl-64-119

**Heifetz Tumor Grasping Forceps**  
**Bayonet Shaped 22 cm - 8 3/4"**  
**Tip Size 3.0 mm Ø**



Jl-64-120

**Heifetz Tumor Grasping Forceps**  
**Bayonet Shaped 22 cm - 8 3/4"**  
**Tip Size 5.0 mm Ø**



Jl-64-121

**Hunt Tumor Grasping Forceps**  
**Angled 20 cm - 8\" Tip Size 5.0**  
**mm Ø**

Tissue Grasping and Tumor Grasping Forceps





Jl-65-101

**Buck Percussion Hammer With  
Needle and Brush 18.5 cm - 7  
1/4" Stainless Steel - Rubber**



Jl-65-102

**Dejerine Percussion Hammer 20  
cm - 8" Stainless Steel - Rubber**



Jl-65-103

**Dejerine Percussion Hammer  
With Needle 21 cm - 8 1/4"  
Stainless Steel - Rubber**



Jl-65-104

**Taylor Mod. USA Percussion  
Hammer 20 cm - 8" Stainless  
Steel - Rubber**



Jl-65-105

**Babinsky Percussion Hammer  
With Needle 22 cm - 8 3/4"  
Stainless Steel - Rubber**



Jl-65-106

**Rabiner Percussion Hammer 23  
cm - 9" Stainless Steel - Rubber**



Jl-65-107

**Wartenberg Sensibility  
Instrument 19 cm - 7 1/2"  
Stainless Steel**



Jl-66-101

**Corona™ Premium ENT  
Diagnostic Set**



Jl-66-102

**Corona™ Premium Otoscope**



Jl-66-103

**Corona™  
Premium Ophthalmoscope**



Jl-66-104

**Apollo™ Standard Universal ENT  
Diagnostic Set**



Jl-66-105

**Apollo™ Standard Otoscope**



Jl-66-106

**Apollo™  
Standard Ophthalmoscope**



Jl-66-107

**Apollo™  
Standard Dermatoscope**



Jl-66-108

**Swisslite LED Illuminate Tongue  
Depressor  
Features**



JI-67-101

**Apollo™ Standard  
Miller Laryngoscope Set With  
Battery Handle**



JI-67-102

**Apollo™ Standard  
Miller Laryngoscope Set With  
Battery Handle**



JI-67-103

**Apollo™ Standard  
Miller Laryngoscope Set With  
Battery Handle**



JI-67-104

**Apollo™ Standard  
Miller Laryngoscope Set With  
Battery Handle**



JI-67-105

**Apollo™ Standard  
Miller Laryngoscope Blade Fig.  
0 - For Babies Working Length  
55 mm**



JI-67-106

**Apollo™ Standard  
Miller Laryngoscope Blade Fig.  
1 - For Children Working Length  
80 mm**



JI-67-107

**Apollo™ Standard  
Miller Laryngoscope Blade Fig.  
2 - For Adolescents Working  
Length 130 mm**



JI-67-108

**Apollo™ Standard  
Miller Laryngoscope Blade Fig.  
4 - For Men Working Length  
180 mm**



JI-67-109

**Apollo™ Standard  
Miller Laryngoscope Blade Fig.  
3 - For Women Working Length  
170 mm**



Jl-67-110

**Apollo™ Standard  
McIntosh-Baby Laryngoscope  
Set With Battery Handle**



Jl-67-111

**Apollo™ Standard McIntosh  
Laryngoscope Set With Battery  
Handle**



Jl-67-112

**Apollo™ Standard McIntosh  
Laryngoscope Set With Battery  
Handle**



Jl-67-113

**Apollo™ Standard McIntosh  
Laryngoscope Set With Battery  
Handle**



Jl-67-114

**Apollo™ Standard McIntosh  
Laryngoscope Blade Fig. 0 - For  
Babies Working Length 55 mm**



Jl-67-115

**Apollo™ Standard  
McIntosh Laryngoscope Blade  
Fig. 1 - For Children Working  
Length 75 mm**



Jl-67-116

**Apollo™ Standard  
McIntosh Laryngoscope Blade  
Fig. 2 - For  
Adolescents Working Length 90  
mm**



Jl-67-117

**Apollo™ Standard  
McIntosh Laryngoscope Blade  
Fig. 3 - For Women Working  
Length 110 mm**



Jl-67-118


**Apollo™ Standard  
McIntosh Laryngoscope Blade  
Fig. 4 - For Men Working Length  
135 mm**






Jl-67-119

**Apollo™ Standard Battery  
Handle**



Jl-67-120

**Apollo™ Standard Spare Bulb  
Large For Blades Fig. 2 to 4**



Jl-67-121

**Apollo™ Standard Spare Bulb  
Small For Blades Fig. 0 to 1**

Standard Laryngoscopes



JI-68-101

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Set  
With Battery Handle**



JI-68-102

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Set  
With Battery Handle**



JI-68-103

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Set  
With Battery Handle**



JI-68-104

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Set  
With Battery Handle**



JI-68-105

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 0 - For Babies Working  
Length 55 mm**



JI-68-106

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 1 - For Children Working  
Length 70 mm**



JI-68-107

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 2 - For Adolescents  
Working Length 80 mm**



JI-68-108

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 3 - For Women Working  
Length 105 mm**



JI-68-109

**Corona™ Premium Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 4 - For Men Working Length  
130 mm**



JI-68-110

**Corona™ Premium Fiber Optic Miller Laryngoscope Set With Battery Handle**



JI-68-111

**Corona™ Premium Fiber Optic Miller Laryngoscope Set With Battery Handle**



JI-68-112

**Corona™ Premium Fiber Optic Miller Laryngoscope Set With Battery Handle**



JI-68-113

**Corona™ Premium Fiber Optic Miller Laryngoscope Set With Battery Handle**



JI-68-114

**Corona™ Premium Fiber Optic Miller Laryngoscope Blade Fig. 0 - For Babies Working Length 55 mm**



JI-68-115

**Corona™ Premium Fiber Optic Miller Laryngoscope Blade Fig. 1 - For Children Working Length 80 mm**



JI-68-116

**Corona™ Premium Fiber Optic Miller Laryngoscope Blade Fig. 2 - For Adolescents Working Length 130 mm**



JI-68-117

**Corona™ Premium Fiber Optic Miller Laryngoscope Blade Fig. 3 - For Women Working Length 170 mm**



JI-68-118

**Corona™ Premium Fiber Optic Miller Laryngoscope Blade Fig. 4 - For Men Working Length 180 mm**



JI-68-119

**Corona™ Premium Fiber  
Optic Battery Handle With 2.5  
Volts Xenon Bulb**



JI-68-120

**Corona™ Premium Fiber Optic  
Battery Handle Slim Profile -  
With 2.5 Volts Xenon Bulb**



JI-68-121

**Corona™ Premium Fiber Optic  
Spare Bulb Xenon**

Fiber Optic Laryngoscopes





JI-69-101

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Set With  
Battery Handle**



JI-69-102

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Set With  
Battery Handle**



JI-69-103

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Set With  
Battery Handle**



JI-69-104

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Set With  
Battery Handle**



JI-69-105

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 0 - For Babies Working  
Length 55 mm**



JI-69-106

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 1 - For Children Working  
Length 70 mm**



JI-69-107

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 2 - For Adolescents Working  
Length 80 mm**



JI-69-108

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 3 - For Women Working  
Length 105 mm**



JI-69-109

**MaxBright™ Fiber Optic  
McIntosh Laryngoscope Blade  
Fig. 4 - For Men Working Length  
130 mm**



JI-69-110

**MaxBright™ Fiber Optic  
Miller Laryngoscope Set With  
Battery Handle**



JI-69-111

**MaxBright™ Fiber Optic  
Miller Laryngoscope Set With  
Battery Handle**



JI-69-112

**MaxBright™ Fiber Optic  
Miller Laryngoscope Set With  
Battery Handle**



JI-69-113

**MaxBright™ Fiber Optic  
Miller Laryngoscope Set With  
Battery Handle**



JI-69-114

**MaxBright™ Fiber Optic  
Miller Laryngoscope Blade Fig. 0  
- For Babies Working Length 55  
mm**



JI-69-115

**MaxBright™ Fiber Optic  
Miller Laryngoscope Blade Fig. 1  
- For Children Working Length 80  
mm**



JI-69-116

**MaxBright™ Fiber Optic  
Miller Laryngoscope Blade Fig. 2  
- For Adolescents Working  
Length 130 mm**



JI-69-117

**MaxBright™ Fiber Optic  
Miller Laryngoscope Blade Fig. 3  
- For Women Working Length  
170 mm**



JI-69-118

**MaxBright™ Fiber Optic  
Miller Laryngoscope Blade Fig. 4  
- For Men Working Length 180  
mm**



JI-69-119

**MaxBright™ Fiber Optic Battery  
Handle With 2.5 Volts LED Bulb**



JI-69-120

**MaxBright™ Fiber Optic Battery  
Handle Slim Profile - With 2.5  
Volts LED Bulb**



JI-69-121

**MaxBright™ Fiber Optic Spare  
Bulb**

LED Fiber Optic Laryngoscopes



Jl-185-101

**Gillies Skin Hook Fig. 1 18 cm - 7"**



Jl-185-102

**Gillies Skin Hook Fig. 2 18 cm - 7"**



Jl-185-103

**Frazier Delicate Hook Sharp 13 cm - 5"**



Jl-185-104

**Frazier Delicate Hook Blunt 13 cm - 5"**



Jl-185-105

**Kleinert-Kutz Delicate Hook Fig. 1 16 cm - 6 1/4"**



Jl-185-106

**Kleinert-Kutz Delicate Hook Fig. 2 16 cm - 6 1/4"**



Jl-185-107

**Kleinert-Kutz Delicate Hook Fig. 3 16 cm - 6 1/4"**



Jl-185-108

**Graefe Skin Hook Sharp 15 cm - 6"**



Jl-185-109

**Graefe Skin Hook Blunt 15 cm - 6"**





Jl-185-110

**Skin Hook Sharp Fig. 1 16 cm -  
6 1/4"**



Jl-185-111

**Skin Hook Sharp Fig. 2 16 cm -  
6 1/4"**



Jl-185-112

**Skin Hook Sharp Fig. 3 16 cm -  
6 1/4"**



Jl-185-113

**Skin Hook Sharp Fig. 4 16 cm -  
6 1/4"**



Jl-185-114

**O-Conner Skin Hook Sharp 16  
cm - 6 1/4"**



Jl-185-115

**Guthrie Skin Retractor  
Sharp 15.5 cm - 6" Width 2.0  
mm**



Jl-185-116

**Guthrie Skin Retractor  
Sharp 15.5 cm - 6" Width 3.0  
mm**



Jl-185-117

**Moldestad Vessel Hook 12 cm -  
4 3/4"**



Jl-185-118

**Brom Vein Hook 19.5 cm - 7  
3/4"**



Jl-185-119

**Brom Vein Hook 19.5 cm - 7  
3/4"**

Delicate Hooks



JI-186-101

**Cushing Retractor / Saddle**  
Hook 20.5 cm - 8" Blade Size 16  
mm



JI-186-102

**Cushing Retractor / Saddle**  
Hook 20.5 cm - 8" Blade Size 14  
mm



JI-186-103

**Cushing Retractor / Saddle**  
Hook 20.5 cm - 8" Blade Size 12  
mm



JI-186-104

**Cushing Retractor / Saddle**  
Hook 20.5 cm - 8" Blade Size 12  
mm



JI-186-105

**Cushing Retractor / Saddle**  
Hook 20.5 cm - 8" Blade Size 10  
mm



JI-186-106

**Cushing Retractor / Saddle**  
Hook 20.5 cm - 8" Blade Size 8  
mm



JI-186-107

**Cushing Retractor / Saddle**  
Hook 20.5 cm - 8"



JI-186-108

**Desmarres Lid Retractor 16 cm**  
- 6 1/4" Blade Width 8.0 mm



JI-186-109

**Desmarres Lid Retractor 16 cm**  
- 6 1/4" Blade Width 10.0 mm



JI-186-110

**Desmarres Lid Retractor 16 cm  
- 6 1/4" Blade Width 12.0 mm**



JI-186-111

**Desmarres Lid Retractor 16 cm  
- 6 1/4" Blade Width 14.0 mm**



JI-186-112

**Desmarres Lid Retractor 16 cm  
- 6 1/4" Blade Width 16.0 mm**

Lid Retractors & Saddle Hooks





Jl-187-101

**Freer Mucosa Hook Sharp 15 cm - 6" Width 3.2 mm**



Jl-187-102

**Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 2.0 mm**



Jl-187-103

**Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 5.0 mm**



Jl-187-104

**Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 7.0 mm**



Jl-187-105

**Joseph Mucosa Hook 16.5 cm - 6 1/2" Width 10.0 mm**



Jl-187-106

**Kilner Alar Hook Sharp 8 cm - 3 1/4"**



Jl-187-107

**Converse Alar Hook Sharp 11.5 cm - 4 1/2"**



Jl-187-108

**Cottle Alar Hook 4 Semi Sharp Prongs 14 cm - 5 1/2"**



Jl-187-109

**Tissue Picker 2 Twisted Sharp Prongs 18.5 cm - 7 1/4" Diameter 1.8 mm Ø**



Jl-188-101

**Joseph Mucus Hook Sharp 15.5 cm - 6"**



Jl-188-102

**Cottle Nasal Hook Sharp 17 cm - 6 3/4"**



Jl-188-103

**Cottle Nasal Hook Sharp - Fig. 1 14.5 cm - 5 3/4"**



Jl-188-104

**Cottle Nasal Hook Sharp - Fig. 2 14.5 cm - 5 3/4"**



Jl-188-105

**Cottle Nasal Hook Sharp - Fig. 3 14.5 cm - 5 3/4"**



Jl-188-106

**Cushing Nerve Hook Blunt - Fig. 1 19 cm - 7 1/2"**



Jl-188-107

**Cushing Nerve Hook Blunt - Fig. 2 19.5 cm - 7 3/4"**



Jl-188-108

**Cushing Nerve Hook Blunt - Fig. 1 28 cm - 11"**



Jl-188-109

**Cushing Nerve Hook Blunt - Fig. 2 28 cm - 11"**



Jl-188-110

**Graefe Nerve Hook Fig. 1 14 cm  
- 5 1/2"**



Jl-188-111

**Graefe Nerve Hook Fig. 2 14 cm  
- 5 1/2"**



Jl-188-112

**Graefe Nerve Hook Fig. 3 14 cm  
- 5 1/2"**

Nerve Hooks



Jl-188-113

**Crile Nerve Hook 20 cm - 8"**



Jl-188-114

**Crile Nerve Hook 20.5 cm - 8"**











Jl-190-101

**Suction Tube 14 cm - 5 1/2"**  
**Diameter 1.5 mm Ø**



Jl-190-102

**Suction Tube 18 cm - 7"**  
**Diameter 1.5 mm Ø**



Jl-190-103

**Suction Tube 14 cm - 5 1/2"**  
**Diameter 2.5 mm Ø**



Jl-190-104

**Frazier Suction Tube 7**  
**Charr. 19.5 cm - 7 3/4"**  
**Diameter 2.3 mm Ø**



Jl-190-105

**Frazier Suction Tube 6**  
**Charr. 19.5 cm - 7 3/4"**  
**Diameter 2.0 mm Ø**



Jl-190-106

**Frazier Suction Tube 8**  
**Charr. 19.5 cm - 7 3/4"**  
**Diameter 2.6 mm Ø**



Jl-190-107

**Fergusson Suction Tube 6**  
**Charr. 18 cm - 7" Diameter 2.0**  
**mm Ø**



Jl-190-108

**Fergusson Suction Tube 8**  
**Charr. 18 cm - 7" Diameter 2.6**  
**mm Ø**



Jl-190-109

**Fergusson Suction Tube 9**  
**Charr. 18 cm - 7" Diameter 3.0**  
**mm Ø**



Jl-190-110

**Wullstein Suction Handle 14 cm - 5 1/2"**



Jl-190-111

**Wullstein Micro Suction Tube 4.5 cm - 1 3/4" Diameter 0.4 mm Ø**



Jl-190-112

**Wullstein Micro Suction Tube 4.5 cm - 1 3/4" Diameter 0.6 mm Ø**



Jl-190-113

**Baron-Schuknecht Suction Tube With Finger Cutt Off Working Length - Diameter 75 mm - 2.0 mm Ø**



Jl-190-114

**House Suction Tube With Finger Cutt Off Suction Diameter - Irrigation Diameter 0.9 mm Ø - 0.7 mm Ø**



Jl-190-115

**House Suction Tube With Finger Cutt Off Suction Diameter - Irrigation Diameter 1.2 mm Ø - 0.9 mm Ø**



Jl-190-116

**Yankauer Suction Tube Fig. 3 13.5 cm - 5 1/4"**



Jl-190-117

**Yankauer Suction Tip**



Jl-190-118

**Poole Suction Tube 22 cm - 8 3/4" Diameter 10.0 mm Ø**





JI-190-119

**Poole Suction Tube 22.5 cm - 8 3/4" Diameter 8.0 mm Ø**



JI-190-120

**Poole-Baby Suction Tube 20 cm - 8" Diameter 5.5 mm Ø**



JI-190-121

**Cooley Suction Tube 36 cm - 14 1/4" Diameter 8.0 mm Ø**



JI-190-122

**Cooley Suction Tube With Perforated Screw Tip 31 cm - 12 1/4" Diameter 10.0 mm Ø**



JI-190-123

**Cooley Suction Tube With Perforated Screw Tip 31 cm - 12 1/4" Diameter 7.0 mm Ø**



JI-190-124

**Cooley Suction Tube With Perforated Screw Tip 36 cm - 14 1/4" Diameter 8.0 mm Ø**



JI-190-125

**DeBakey Suction Tube 27.5 cm - 10 3/4" Diameter 9.0 mm Ø**



JI-190-126

**Finsterer Suction Tube 32 cm - 12 3/4" 15.0 mm Ø**



JI-190-127

**Nuboer Suction Tube Round Opening 21 cm - 8 1/4" Diameter 6.0 mm Ø**



JI-190-128

**Byrd Suction Tube 21.5 cm - 8  
1/2" Diameter 3.5 mm Ø**



JI-190-129

**Cogswell Suction Tube With  
Finger Cutt Off 24.5 cm - 9 3/4"  
Diameter 2.0 mm Ø**

Suction Tubes



JL-191-101

**Universal Trocar With 4 Trocar  
Tips Diameter 1.7 mm, 3.0 mm,  
4.5 mm and 6.0 mm 14.5 cm - 5  
3/4"**



JL-191-102

**Nelson Trocar 22 FG - Fig. 1 20  
cm - 8" Diameter 7.5 mm Ø**



JL-191-103

**Nelson Trocar 33 FG - Fig. 3 20  
cm - 8" Diameter 11.0 mm Ø**

Trocars



JL-191-104

**Nelson Trocar 28 FG - Fig. 2 20  
cm - 8" Diameter 9.5 mm Ø**



JL-191-105

**Ochsner Trocar 14 cm - 5 1/2"  
Diameter 4.0 mm Ø**



JL-191-106

**Ochsner Trocar 14 cm - 5 1/2"  
Diameter 4.6 mm Ø**



JL-191-107

**Ochsner Trocar 14 cm - 5 1/2"  
Diameter 5.3 mm Ø**



JL-191-108

**Ochsner Trocar 14 cm - 5 1/2"  
Diameter 6.0 mm Ø**



JL-191-109

**Ochsner Trocar 14 cm - 5 1/2"  
Diameter 6.6 mm Ø**



Jl-191-110

**Ochsner Trocar 14 cm - 5 1/2"**  
**Diameter 8.0 mm Ø**



Jl-191-111

**Ochsner Trocar 14 cm - 5 1/2"**  
**Diameter 7.3 mm Ø**



Jl-191-112

**Ochsner Trocar 14 cm - 5 1/2"**  
**Diameter 8.6 mm Ø**

Trocars





JI-192-101

**Rollet Fine Wound Retractor 2  
Sharp Prongs 13 cm - 5" Width  
4.5 mm**



JI-192-102

**Rollet Fine Wound Retractor 4  
Sharp Prongs 13.5 cm - 5 1/4"  
Width 7.0 mm**



JI-192-103

**Rollet Fine Wound Retractor 4  
Blunt Prongs 13.5 cm - 5 1/4"  
Width 7.0 mm**



JI-192-104

**Knapp Fine Wound Retractor 4  
Sharp Prongs 13.5 cm - 5 1/4"**



JI-192-105

**Knapp Fine Wound Retractor 4  
Blunt Prongs 13.5 cm - 5 1/4"**



JI-192-106

**McIndoe Wound Retractor 2  
Sharp Prongs 15 cm - 6" Width  
9.0 mm**



JI-192-107

**McIndoe Wound Retractor 2  
Sharp Prongs 15 cm - 6" Width  
10.0 mm**



JI-192-108

**Wound Retractor 1 Sharp  
Prong 16.5 cm - 6 1/2"**



JI-192-109

**Wound Retractor 2 Sharp  
Prongs - Small Curve 16.5 cm -  
6 1/2" Width 4.2 mm**



JI-192-110

**Wound Retractor 3 Sharp  
Prongs - Small Curve 16.5 cm -  
6 1/2" Width 7.0 mm**



JI-192-111

**Wound Retractor 4 Sharp  
Prongs - Small Curve 16.5 cm -  
6 1/2" Width 9.5 mm**



JI-192-112

**Wound Retractor 1 Blunt  
Prong 16.5 cm - 6 1/2"**



JI-192-113

**Wound Retractor 2 Blunt Prongs  
- Small Curve 16.5 cm - 6 1/2"  
Width 4.2 mm**



JI-192-114

**Wound Retractor 3 Blunt Prongs  
- Small Curve 16.5 cm - 6 1/2"  
Width 7.0 mm**



JI-192-115

**Wound Retractor 4 Blunt Prongs  
- Small Curve 16.5 cm - 6 1/2"  
Width 9.5 mm**



JI-192-116

**Wound Retractor Flexible - 2  
Sharp Prongs 16 cm - 6 1/4"**



JI-192-117

**Wound Retractor Flexible - 3  
Sharp Prongs 16 cm - 6 1/4"**



JI-192-118

**Wound Retractor Flexible - 4  
Sharp Prongs 16 cm - 6 1/4"**



JI-192-119

**Wound Retractor Flexible - 2  
Blunt Prongs 16 cm - 6 1/4"**



JI-192-120

**Wound Retractor Flexible - 3  
Blunt Prongs 16 cm - 6 1/4"**



JI-192-121

**Wound Retractor Flexible - 4  
Blunt Prongs 16 cm - 6 1/4"**

Wound Retractors



JI-192-122

**Wound Retractor 4 Sharp  
Prongs 16.5 cm - 6 1/2"**



JI-192-123

**Wound Retractor 4 Blunt  
Prongs 16.5 cm - 6 1/2"**



JI-192-124

**Wound Retractor 1 Sharp  
Prong 16.5 cm - 6 1/2"**



Jl-193-101

**Senn-Miller Retractor Blunt 16 cm - 6 1/4" Blade Size 20 x 7 mm**



Jl-193-102

**Senn-Miller Retractor Sharp 16 cm - 6 1/4" Blade Size 20 x 7 mm**



Jl-193-103

**Ragnell-Davis Retractor 14 cm - 5 1/2" Blade Size 8 x 4 mm - 15 x 5.5 mm**



Jl-193-104

**Klapp Retractor 2 Blunt Prongs 17.5 cm - 7" Blade Size 13 x 6 mm**



Jl-193-105

**Klapp Retractor 3 Blunt Prongs 17.5 cm - 7" Blade Size 12 x 11 mm**



Jl-193-106

**Klapp Retractor 4 Blunt Prongs 17.5 cm - 7" Blade Size 15 x 11 mm**



Jl-193-107

**Cushing-Kocher Retractor 17.5 cm - 7" Blade Size 10 x 10 mm**



Jl-193-108

**Retractor 4 Sharp Prongs 16.5 cm - 6 1/2"**



Jl-193-109

**Koenig Retractor 18 cm - 7" Blade Size 12 x 16 mm**





JI-193-110

**Retractor 4 Blunt Prongs 16.5 cm - 6 1/2"**



JI-193-111

**Koenig Retractor 18 cm - 7" Blade Size 13 x 13 mm**



JI-193-112

**Volkmann Retractor 1 Sharp Prong - Hollow Handle 21.5 cm - 8 1/2"**



JI-193-113

**Zenker Retractor 30.5 cm - 12" Blade Size 100 x 25 mm**



JI-193-114

**Volkmann Retractor 2 Sharp Prongs - Hollow Handle 21.5 cm - 8 1/2"**



JI-193-115

**Zenker Retractor 28.5 cm - 11 1/4" Blade Size 85 x 20 mm**



JI-193-116

**Volkmann Retractor 3 Sharp Prongs - Hollow Handle 21.5 cm - 8 1/2"**



JI-193-117

**Zenker Retractor 27.5 cm - 10 3/4" Blade Size 75 x 12 mm**



JI-193-118

**Volkmann Retractor 4 Sharp Prongs - Hollow Handle 21.5 cm - 8 1/2"**



JI-193-119

**Zenker Retractor 26.5 cm - 10 1/2" Blade Size 60 x 12 mm**



JI-193-120

**Volkmann Retractor 6 Sharp Prongs - Hollow Handle 21.5 cm - 8 1/2"**



JI-193-121

**Zenker Retractor 25 cm - 9 3/4" Blade Size 35 x 8 mm**



JI-193-122

**Kirsch Retractor Right 21 cm - 8 1/4" Blade Size 60 x 11 mm**



JI-193-123

**Kirsch Retractor Left 21 cm - 8 1/4" Blade Size 60 x 11 mm**



JI-193-124

**Volkmann Retractor 1 Blunt Prong - Hollow Handle 21.5 cm - 8 1/2"**



JI-193-125

**Wound Retractor 2 Blunt Prongs - Large Curve 16.5 cm - 6 1/2" Width 7.5 mm**



JI-193-126

**Volkmann Retractor 2 Blunt Prongs - Hollow Handle 21.5 cm - 8 1/2"**



JI-193-127

**Volkmann Retractor 3 Blunt Prongs - Hollow Handle 21.5 cm - 8 1/2"**



JI-193-128

**Volkmann Retractor 4 Blunt  
Prongs - Hollow Handle 21.5 cm  
- 8 1/2"**



JI-193-129

**Volkmann Retractor 6 Blunt  
Prongs - Hollow Handle 21.5 cm  
- 8 1/2"**



JI-193-130

**Volkmann Retractor 1 Sharp  
Prong 22 cm - 8 3/4"**



JI-193-131

**Volkmann Retractor 2 Sharp  
Prongs 22 cm - 8 3/4"**



JI-193-132

**Volkmann Retractor 3 Sharp  
Prongs 22 cm - 8 3/4"**



JI-193-133

**Volkmann Retractor 4 Sharp  
Prongs 22 cm - 8 3/4"**



JI-193-134

**Volkmann Retractor 6 Sharp  
Prongs 22 cm - 8 3/4"**



JI-193-135

**Volkmann Retractor 1 Blunt  
Prong 22 cm - 8 3/4"**



JI-193-136

**Volkmann Retractor 2 Blunt  
Prongs 22 cm - 8 3/4"**



Jl-193-137

**Volkman Retractor 3 Blunt**  
Prongs 22 cm - 8 3/4"



Jl-193-138

**Volkman Retractor 4 Blunt**  
Prongs 22 cm - 8 3/4"



Jl-193-139

**Volkman Retractor 6 Blunt**  
Prongs 22 cm - 8 3/4"



Jl-193-140

**Volkman Retractor 1 Blunt**  
Prong 21.5 cm - 8 1/2"



Jl-193-141

**Kocher-  
Langenbeck Retractor 21.5 cm -  
8 1/2" Blade Size 35 x 8 mm**



Jl-193-142

**Kocher-  
Langenbeck Retractor 21.5 cm -  
8 1/2" Blade Size 25 x 6 mm**



Jl-193-143

**Volkman Retractor 1 Sharp**  
Prong 21.5 cm - 8 1/2"



Jl-193-144

**Kocher Retractor 23 cm - 9"**  
Blade Size 75 x 40 mm



Jl-193-145

**Kocher Retractor 1 Sharp**  
Prong 22 cm - 8 3/4"





JI-193-146

**Kocher Retractor 2 Sharp  
Prongs 22 cm - 8 3/4"**



JI-193-147

**Kocher Retractor 23 cm - 9"  
Blade Size 75 x 30 mm**



JI-193-148

**Kocher Retractor 23 cm - 9"  
Blade Size 61 x 25 mm**



JI-193-149

**Kocher Retractor 3 Sharp  
Prongs 22 cm - 8 3/4"**



JI-193-150

**Kocher Retractor 23 cm - 9"  
Blade Size 61 x 20 mm**




JI-193-151

**Kocher Retractor 23 cm - 9"  
Blade Size 41 x 18 mm**



JI-193-152

**Israel Retractor 7 Blunt  
Prongs 25.5 cm - 10" Blade Size  
70 x 70 mm**




JI-193-153

**Ollier Retractor 3 Blunt  
Prongs 23 cm - 9" Blade Size 36  
x 30 mm**



JI-193-154

**Israel Retractor 6 Blunt  
Prongs 25.5 cm - 10" Blade Size  
50 x 60 mm**



Jl-193-155

**Ollier Retractor 4 Blunt**  
Prongs 23 cm - 9" Blade Size 36  
x 60 mm




Jl-193-156

**Israel Retractor 5 Blunt**  
Prongs 25.5 cm - 10" Blade Size  
45 x 50 mm




Jl-193-157

**Wassmund Retractor 20.5 cm -**  
8" Blade Size 33 x 21 mm




Jl-193-158

**Israel Retractor 4 Blunt**  
Prongs 25.5 cm - 10" Blade Size  
40 x 40 mm



Jl-193-159

**Koerte Retractor 8 Blunt**  
Prongs 24.5 cm - 9 3/4" Blade  
Size 40 x 30 mm



Jl-193-160

**Koerte Retractor 8 Sharp**  
Prongs 24.5 cm - 9 3/4" Blade  
Size 40 x 30 mm




Jl-193-161

**Kocher Retractor 21.5 cm - 8**  
1/2" Blade Size 22 x 12 mm



Jl-193-162

**Kocher Retractor 21.5 cm - 8**  
1/2" Blade Size 24 x 14 mm




Jl-193-163

**Hoesel Retractor 26 cm - 10**  
1/4" Blade Size 83 x 30 mm



Jl-193-164

**Deaver Retractor Fig. 6 36 cm -  
14 1/4" Blade Width 25 mm**




Jl-193-165

**Brunner Retractor 25.5 cm - 10"  
Blade Size 104 x 25 mm**



Jl-193-166

**Deaver Retractor Fig. 5 33 cm -  
13" Blade Width 25 mm**




Jl-193-167

**Brunner Retractor 25.5 cm - 10"  
Blade Size 124 x 30 mm**




Jl-193-168

**Deaver Retractor Fig. 7 30 cm -  
11 3/4" Blade Width 38 mm**




Jl-193-169

**Brunner Retractor 25.5 cm - 10"  
Blade Size 130 x 30 mm**



Jl-193-170

**Brunner Retractor 25.5 cm - 10"  
Blade Size 144 x 30 mm**




Jl-193-171

**Brunner Retractor 25.5 cm - 10"  
Blade Size 160 x 30 mm**



Jl-193-172

**Deaver Retractor Fig. 8 31 cm -  
12 1/4" Blade Width 50 mm**



JI-193-173

**Brunner Retractor 25.5 cm - 10"**  
**Blade Size 164 x 30 mm**



JI-193-174

**Deaver Retractor Fig. 9 31.5 cm**  
**- 12 1/2" Blade Width 75 mm**



JI-193-175

**Brunner Retractor 25.5 cm - 10"**  
**Blade Size 82 x 20 mm**



JI-193-176

**Kelly Retractor 26 cm - 10 1/4"**  
**Blade Size 62 x 50 mm**



JI-193-177

**Kelly Retractor 26 cm - 10 1/4"**  
**Blade Size 70 x 65 mm**



JI-193-178

**Kelly Retractor Stainless**  
**Steel 26 cm - 10 1/4" Blade**  
**Size 90 x 76 mm**



JI-193-179

**Tuffier Retractor 21.5 cm - 8**  
**1/2" Blade Size 58 x 64 mm**



JI-193-180

**Morris Retractor 24.5 cm - 9**  
**3/4" Blade Size 70 x 40 mm**



JI-193-181

**Morris Retractor 24.5 cm - 9**  
**3/4" Blade Size 70 x 50 mm**





Jl-193-182

**Morris Retractor 24.5 cm - 9 3/4" Blade Size 70 x 65 mm**



Jl-193-183

**Koerte Retractor 30 cm - 11 3/4" Blade Width 46 mm**



Jl-193-184

**Cooley Retractor Right - Serrated 25 cm - 9 3/4" Blade Size 48 x 48 mm**



Jl-193-185

**Mikulicz Liver Retractor 25 cm - 9 3/4" Blade Size 90 x 35 mm**



Jl-193-186

**Mikulicz Liver Retractor 25 cm - 9 3/4" Blade Size 125 x 50 mm**



Jl-193-187

**Mikulicz Liver Retractor 25 cm - 9 3/4" Blade Size 160 x 50 mm**



Jl-193-188

**Doyen Retractor 28 cm - 11" Blade Size 48 x 120 mm**



Jl-193-189

**Doyen Retractor 28 cm - 11" Blade Size 48 x 135 mm**



Jl-193-190

**Chain with Weight 79.0 cm - 31"**



Jl-193-191

**Kocher Retractor 24 cm - 9 1/2"**  
**Blade Size 73 x 25 mm**



Jl-193-192

**Kocher Retractor 25 cm - 9 3/4"**  
**Blade Size 80 x 45 mm**



Jl-193-193

**Kocher Retractor 25 cm - 9 3/4"**  
**Blade Size 80 x 55 mm**



Jl-193-194

**Kocher Retractor 25 cm - 9 3/4"**  
**Blade Size 80 x 65 mm**



Jl-193-195

**Doyen-Stillé Retractor 26 cm - 10 1/4"**  
**Blade Size 87 x 28 mm**



Jl-193-196

**Doyen-Stillé Retractor 26 cm - 10 1/4"**  
**Blade Size 97 x 31 mm**



Jl-193-197

**Doyen-Stillé Retractor 26 cm - 10 1/4"**  
**Blade Size 117 x 52 mm**



Jl-193-198

**Fritsch Retractor 23.5 cm - 9 1/4"**  
**Blade Size 37 x 31 mm**



Jl-193-199

**Fritsch Retractor 24 cm - 9 1/2"**  
**Blade Size 38 x 40 mm**



Jl-193-200

**Fritsch Retractor 25 cm - 9 3/4"**  
**Blade Size 43 x 51 mm**



Jl-193-201

**Fritsch Retractor 25 cm - 9 3/4"**  
**Blade Size 49 x 61 mm**



Jl-193-202

**Fritsch Retractor 25.5 cm - 10"**  
**Blade Size 48 x 75 mm**



Jl-193-203

**Fritsch Retractor 26.5 cm - 10 1/2"**  
**Blade Size 65 x 85 mm**



Jl-193-204

**Collin-Hartmann Retractor 15 cm - 6"**



Jl-193-205

**Sternberg Retractor 12 cm - 4 3/4"**



Jl-193-206

**Sternberg Retractor 13 cm - 5"**



Jl-193-207

**Sternberg Retractor 13.5 cm - 5 1/4"**



Jl-193-208

**Roux-Mini Retractor 14.5 cm - 5 3/4"**  
**Blade Size 6 x 7 mm - 7 x 8 mm**



193-209

**Baby-Roux Retractor 12.5 cm -  
5" Blade Size 17 x 17 mm - 24 x  
21 mm**



193-210

**Roux Retractor Set of Fig. 1,  
Fig. 2 and Fig. 3**



193-211

**Roux Retractor Fig. 1 14 cm - 5  
1/2"**

General Retractors



# CONTACT US

## Address:

PATESIR JAMMU ROAD P.O.BOX DALLOWALI,  
Sialkot 51310 Pakistan

**Phone:** +92-52-3206311

**FAX:** +92-52-3205458

**Cell:** +92-331-6116799 

**Skype:** julespak.int

**E-mail:** [info@julespak.com](mailto:info@julespak.com)

**Web:** [www.julespak.com](http://www.julespak.com)